STORZ

KARL STORZ — ENDOSKOPE

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INSTRUCTION MANUAL
UP 501 S-PILOT™ Tube Squeeze Valve

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**1 Important information for
users of KARL STORZ
devices**

**It is recommended to check the suitability of
the product for the intended procedure prior
to use.**

Thank you for your expression of confidence in
the KARL STORZ brand name. Like all of our
products, this product too is the result of years
of experience and careful manufacturing. You
and your organization have decided in favor of a
modern high quality product from KARL STORZ.

This instruction manual is intended to serve as
an aid in the proper installation, connection and
operation of the tube squeeze valve S-PILOT™
UP 501. All of the necessary details and actions
are clearly explained.

Please read these instructions carefully. Keep this
manual in a convenient and conspicuous location
and in its protective casing close to the device.

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CAUTION: Federal (USA) law restricts this device
to sale by or on the order of a physician.

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**3 Controls, displays,
connectors, and their uses**

- ① On/Off switch with status LED (orange or green)
- ② Button for opening/closing the tube squeeze valve with status LED
'green' = open
'off' = closed
- ③ Tube squeeze valve
- ④ Connection socket for footswitch 20 0141 30 or control cable 27 6100 70 or 20 7010 70 for connecting the UNIDRIVE® S III or CALCUSON
- ⑤ Connection socket for 3.5 mm connecting cable UP 004 or UP 005 (for optional remote control of an HF surgical unit)
- ⑥ Slot for SD memory card (for software updates)
- ⑦ SCB connectors
- ⑧ Power cord socket
- ⑨ Potential equalization connector

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	Follow instructions for use
	ON/OFF (Standby)
	Tube squeeze valve button
	Applied part of type CF
	Connection socket for footswitch
	Potential equalization connector
	Remote control
	Alternating current
	Electronic information product pollution control (China RoHS)
	This device has been labelled in Accordance with the European Directive on waste electrical and electronic equipment (WEEE).
	Manufacturer

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
For the meanings of the symbols printed on the label or packaging, please refer to the 'Packaging symbols' accompanying instruction leaflet, mat. no. 96216316 DF. This can be downloaded from www.karlstorz.com.

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
	Manufacturer
	Catalogue number
	Do not reuse
	Do not resterilize
	CE mark
	Standard compliance According to IEC 60601-1: Degree of protection against electric shocks – Applied part type CF
	Do not use if package is damaged
	Consult instructions for use
	In the USA, federal law restricts this device to sale or use by or on the order of a physician
	Store in a dry place
	Keep away from sunlight
	Sterilized with ethylene oxide
	Use-by date
	Batch code

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1	Important information for users of KARL STORZ devices	III
2	Images of the equipment	IV
3	Controls, displays, connectors, and their uses	V
4	Symbols employed	VI
4.1	Symbols on the device	VI
4.2	Symbols on label and packaging of the device	VI
4.3	Tubing set label symbols	VII
5	General information	3
5.1	Description of the device	3
5.2	Property rights	5
5.3	Free RTOS Licensing	6
6	Safety instructions	7
6.1	Explanation of warnings and cautions	7
6.2	Intended use	11
6.2.1	Indications	11
6.2.2	Contraindications	12
6.3	User qualification	12
6.4	User profile of physician and assistants	12
6.5	Intended conditions of use	13
6.5.1	Use	13
6.5.2	Other intended conditions	13
6.6	Safety precautions at the installation site	14
6.7	Safety precautions when operating the device	15
6.8	Safety features	15
6.8.1	Self-test	15
6.8.2	Motor current monitoring	15
7	Installation and operating instructions	16
7.1	Unpacking the equipment	16
7.2	Basic equipment	16
7.3	Installing and connecting the device	16
7.3.1	Connecting the ground line	17
7.3.2	Connecting the power cord	17
7.3.3	Connecting the footswitch for direct control	18
7.3.4	Combining the HF device with the S-PILOT™ in the OR1™ system via the KARL STORZ-SCB – automatic smoke evacuation	18
7.3.5	Combining the HF device with the S-PILOT™ via the KARL STORZ-SCB – automatic smoke evacuation as a standalone solution	18

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7.3.6	Combining the HF device with the S-PILOT™ via the connecting cable UP 004/UP 005 – automatic smoke evacuation as a standalone solution	20
7.3.7	Combining the UNIDRIVE® S III with the S-PILOT™ – suction of fluids	20
7.3.8	Combining the CALCUSON with the S-PILOT™ – suction of fluids	20
7.4	Commissioning	21
7.4.1	Inserting the suction tube set	21
7.4.2	Test for proper functioning	23
7.4.3	Decommissioning	24
8	Maintenance	25
8.1	Reprocessing	25
8.1.1	Reprocessing of the UP 501	26
8.1.2	Reprocessing of suction tube set for single use, sterile	26
8.2	Maintenance and safety check	27
8.2.1	Maintenance	27
8.2.2	Safety checks/ repeat inspections as per IEC 62353	27
8.3	Servicing and repair	28
8.4	Disposal	28
8.5	Repair program	29
8.6	Important information	29
8.7	Limitation of liability	30
8.8	Warranty	30
9	Technical description	31
9.1	Acoustic signals	31
9.2	Troubleshooting	32
9.3	Technical data	34
9.3.1	Standard compliance	35
9.3.2	Directive compliance	35
9.4	Technical documentation	36
10	Spare parts, recommended accessories	37
10.1	Spare parts/accessories	37
10.2	Accessories	37
11	Electromagnetic Compatibility (EMC) information	38
12	Subsidiaries	54

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5 General information

5.1 Description of the device

The tube squeeze valve S-PILOT™ UP 501 with integrated SCB module is used to clamp off or release a suction tube connected to a vacuum reservoir. The underpressure is generated by a suction pump or via the central vacuum supply of the hospital.

The tube squeeze valve can be operated either directly or indirectly.

When operated directly, the footswitch 200141 30 is connected to the footswitch socket on the S-PILOT™. When the footswitch is actuated, the valve opens for the duration of actuation. The valve closes when the activation signal stops.

The tube squeeze valve can also be opened and closed using the button ② on the S-PILOT™ (e.g., by unsterile OR staff).

When operated indirectly, the tube squeeze valve is activated via a device connected to the S-PILOT™. The devices are connected to each other either via the S-PILOT™ footswitch socket, the KARL STORZ-SCB interface or via a connecting cable (UP 004/UP 005).

a) Automatic smoke evacuation using HF surgical unit

The following HF surgical units are supported:

- KARL STORZ AUTOCON® III 400
- KARL STORZ AUTOCON® II 400
- ERBE Vio 300D
- COVIDIEN Force Triad
- COVIDIEN Valleylab Force FX
- KLS Martin maxium ME402
- ConMed System 5000
- ConMed System 2450

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* Please note the following exception: If the AUTOCON® III 400 is integrated in the OR1™ system, the connecting cable UP 004 is additionally required.

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With automatic support, a difference is made between two cases:

- The supported HF device which is used is integrated in the OR1™ system. It is then sufficient for the connection S-PILOT™ – HF device that the S-PILOT™ is integrated into the OR1™ system (via SCB).*
- The S-PILOT™ and the supported HF device (+ accessories) are to be used as an independent solution. In this case, observe the following table.

Device	Accessories
KARL STORZ AUTOCON® III 400	Connecting cable UP 004
KARL STORZ AUTOCON® II 400	SCB connecting cable
ERBE Vio 300D	Set for SCB/ERBE Vio Interface 20 0900 13
COVIDIEN Force Triad	Connecting cable UP 005
COVIDIEN Valleylab Force FX	Connecting cable UP 005
ConMed System 5000	Connecting cable UP 004
ConMed System 2450	Connecting cable UP 004
KLS Martin maxium ME402	Set for SCB/KLS Martin maxium 20 0912 90**
**the KARL STORZ-SCB Interface Control Set is also additionally required	

Smoke gas is produced when operating the surgical unit. If the S-PILOT™ and the HF surgical unit are connected via the SCB interface, the tube squeeze valve on the S-PILOT™ is automatically opened when the HF surgical unit is activated so that smoke gas is suctioned away and the surgeon has an excellent view of the operation site. As soon as the HF signal on the tube squeeze valve stops, the valve closes again after a time delay of 3 seconds.

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b) Suctioning fluid from the shaver attachment using the motor system UNIDRIVE® S III

The signal of a two-stage single-pedal footswitch connected to the UNIDRIVE® S III is looped through the control cable 27 6100 70 to the S-PILOT™. If the surgeon activates the first stage of the footswitch, the valve on the S-PILOT™ opens so that fluid can be suctioned via the vacuum. If the footswitch is pressed into the second stage, the motor system is additionally activated. As the valve also stays open in the second footswitch stage, tissue fragments which have been removed from the shaver, including the rinsing liquid can be suctioned away.

i **NOTE: The KARL STORZ S-PILOT™ together with a KARL STORZ urologic morcellator or hysteroscopic shaver is not authorized for users inside the U.S.**

c) Suctioning fluid when using with the CALCUSON ultrasonic generator

The signal of a two-stage single-pedal footswitch connected to the CALCUSON is looped through the control cable 20 7010 70 to the S-PILOT™. If the surgeon activates the first stage of the footswitch, the valve on the S-PILOT™ opens so that fluid can be suctioned via the vacuum. If the footswitch is pressed into the second stage, the CALCUSON is additionally activated. As the valve also stays open in the second footswitch stage, the stone fragments including the fluid can be suctioned away.

5.2 Property rights

This product is protected in the USA by (at least one of the following) US Patent(s) 5,788,688; 6,397,286; 6,484,221; 6,824,539.

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5.3 Free RTOS Licensing

i **NOTE:** This product uses the FreeRTOS.org real time kernel. The FreeRTOS.org source code can be obtained by visiting <http://www.FreeRTOS.org>.

The FreeRTOS source code is licensed by the GNU General Public License (GPL) with an exception.

The full text of the GPL is available here: <http://www.freertos.org/license.txt>.

The text of the exception is available on FreeRTOS official website: <http://www.FreeRTOS.org> – License and Warranty Page

The exceptions permits the source code of applications that use FreeRTOS solely through the API published on this WEB site to remain closed source, thus permitting the use of FreeRTOS in commercial applications without necessitating that the whole application be open sourced.

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6 Safety instructions

Safety instructions are measures intended to protect the user and patients from the risks which could arise through the use of the system.

6.1 Explanation of warnings and cautions

Please read this instruction manual through carefully and follow its instructions exactly. The words **Warning**, **Caution**, and **Note** convey special meanings. Wherever they are used in this instruction manual, the following text must be carefully reviewed to ensure the safe and effective operation of this device. In addition, to make the signal words stand out more clearly, they are accompanied by a pictogram.

WARNING: A Warning indicates that the personal safety of the patient or physician may be endangered. Failure to observe a Warning could result in injury to the patient or physician.

CAUTION: A Caution indicates that particular service procedures or precautions must be followed to avoid possible damage to the device.

NOTE: A Note provides special information regarding the operation of the product, or clarifies important information.

WARNING: Read this instruction manual thoroughly and be familiar with its contents prior to using this device. In particular, read the chapter on safety instructions to avoid putting your patients, your personnel, or yourself at risk.

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Safety instructions
Warnings and cautions

-  **WARNING:** The electrical installation of the operating room in which the device is connected and operated must comply with the applicable IEC standards.
-  **WARNING:** Set up the device out of reach of patients.
-  **WARNING:** The instruction manuals and interface specifications for medical devices and/or system components used in combination must be observed precisely.
-  **WARNING:** Combinations of medical devices are only assured to be safe if
 - they are identified as such in the respective instruction manuals or
 - the intended use and interface specifications of the devices used in combination permit this (cf. IEC 60601-1-1 or para 16 of the 3rd edition of the IEC 60601-1).
-  **WARNING:** Grounding reliability can only be achieved when the device is connected to a 'Hospital only' or 'Hospital Grade' outlet. Check the plugs and cables regularly and do not use if damaged.
-  **WARNING:** To avoid the risk of electric shock, this equipment must only be connected to a supply mains with protective earth.
-  **WARNING:** No responsibility for the safe operation of the device can be accepted if any tubing systems other than those specified are used.
-  **WARNING:** Before you start the intervention, make sure that the tubing set is not leaking or kinked and that the ends of the tubes are properly secured.
-  **WARNING:** Always unplug the device from the power supply before carrying out any cleaning and maintenance work.

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Safety instructions
Warnings and cautions

-  **WARNING:** Danger of electric shock! Do not open the device. Arrange for service work to be carried out by the manufacturer or by personnel authorized by the manufacturer only. Any opening of the device by unauthorized persons will void warranty claims.
-  **WARNING:** The signal input and signal output parts of this device are designated by the manufacturer for exclusive connection to equipment which complies with IEC 60601-1.
-  **WARNING:** Only operate the device under the specified environmental conditions.
-  **WARNING:** The device may be subject to sudden failures during use. Keep a spare device at hand or alternative procedures in mind if the intervention requires continuously controllable suction.
-  **WARNING:** For safety reasons, do not simultaneously touch the device output sockets and the patient.
-  **WARNING:** Risk of infection: Incorrectly reprocessed medical devices expose patients, users and third parties to a risk of infection as well as the risk that the medical device may malfunction. Observe the 'Cleaning, Disinfection, Care, and Sterilization of KARL STORZ Instruments' instructions and the accompanying documentation.
-  **WARNING:** Check the packaging of the tube set for damage! If the packaging is damaged, the sterility of the tube set may be impaired and the tube set must not be used.
-  **WARNING:** Follow the country-specific regulations/laws for the disposal of accessories.

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WARNING: Avoid allowing fluids to enter the housing. Do not store liquids on or directly above the device. If, despite the precautions, liquid enters the device, sufficient time should be planned for evaporation (this also applies for the formation of condensation).

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Before using the device on the patient for the first time it is imperative that you be acquainted with how the device operates and is controlled.

6.2 Intended use

Valves are used to regulate the volume flow during diagnostic or therapeutic interventions.

6.2.1 Indications

Valves with corresponding accessories are suitable for all diagnostic and therapeutic interventions which require controlled fluid transportation for successful operation.

Use of the device in fields other than those indicated above is not allowed for safety reasons. The S-PILOT™ may only be operated with tubing sets and accessories which have been designated as suitable for the device by KARL STORZ.

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WARNING: No modification of this equipment is allowed.

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WARNING: Additional devices which are connected to electrical medical devices must comply with the relevant IEC or ISO standards (e.g., IEC 60950 for data processing devices). Furthermore, all configurations must comply with the requirements for medical electrical systems (see IEC 60601-1-1 or para 16 of the 3rd edition of the IEC 60601-1). Anybody connecting additional devices to electrical medical devices is configuring a medical system, and is therefore responsible for ensuring that the system complies with the standard requirements for medical electrical systems. It should be noted that the local legal regulations have priority over the above-mentioned requirements. If in doubt, contact your representative on site or the technical customer service.

6.2.2 Contraindications

No contraindications relating directly to the use of this product are currently known.

The use of this device is contraindicated if, in the opinion of an experienced physician, the patient's health is endangered through its use.

6.3 User qualification

The S-PILOT™ may only be used by physicians and medical assistants who have a corresponding specialized qualification and who have been instructed in use of the device.

6.4 User profile of physician and assistants

Recognized medical skills in the field of application (specialist physician, qualified medical staff).

Adequate powers of comprehension to rationally assess the situation in hand.

Familiarity with the use of technical apparatus.

Suitable command of one of the languages used by the device and in the instruction manual.

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Must have been trained thoroughly at least once in the use of the device.

No physical impairments which limit or prevent perception of acoustic or visual alarm signals.

6.5 Intended conditions of use

6.5.1 Use

The device is intended for use in hospitals and doctors' offices. The technical data and ambient conditions are described in the instruction manual.

6.5.2 Other intended conditions

- Frequency of use: one or more times a day
- Length of use: from a few minutes to several hours a day.
- Place of installation: positioning on a level, vibration-free surface
- Mobility: can be moved if positioned on a cart
- Combination: can be used on the patient at the same time as other equipment required for the operations
- Control: can be controlled via the KARL STORZ-SCB.
- The device is placed outside the sterile area.

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6.6 Safety precautions at the installation site

The device may only be used in medical rooms whose electrical systems have been installed in accordance with applicable national regulations.

! WARNING: **DO NOT** use the device in the presence of flammable anesthetics.

! WARNING: The device is not intended for use in potentially explosive zones. Do not operate the device within the demarcated hazard zone while explosive anesthetic gases are in use.

! WARNING: The device must not be used in oxygenated environments.

This also applies for easily combustible and explosive chemicals, e.g. skin disinfectants and fast-acting surface disinfectants.

! WARNING: The device is equipped with a connector for attaching a ground line. It should be connected up in accordance with the applicable national regulations.

! WARNING: The use of additional multiple socket outlets or an extension cable with the system is not permitted.

! WARNING: The use of freely accessible multiple socket outlets is not permitted.

! WARNING: Multiple socket outlets must not be placed on the floor.

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6.7 Safety precautions when operating the device

It is the user's responsibility to make sure the equipment is safe and operates properly before use.

This device has been tested and found to comply with the EMC limits of the Medical Device Directive 93/42/EEC (EN 55011 Class B and IEC 60601-1-2).

Please observe the Electromagnetic Compatibility information in section 11.

6.8 Safety features

The S-PILOT™ has the following safety features:

6.8.1 Self-test

A self-check procedure is performed each time the device is switched on. If the test is negative, this is indicated by the blinking of the status LED and an acoustic alarm.

Its intended use is no longer possible.

The following components have been tested:

- Functionality of the tube squeeze valve

6.8.2 Motor current monitoring

The control of the motor drive of the S-PILOT™ tube squeeze valve contains a temperature and current monitor which switches off when the motor is overheated or overloaded.

i NOTE: A temporary error may be rectified by pressing the button ② on the front panel.

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**7 Installation and operating
instructions**

7.1 Unpacking the equipment

Carefully remove the S-PILOT™ and the accessories from their packaging. Check the delivery for missing items and possible damage.

Please file any complaints with the manufacturer or supplier immediately.

If possible, retain the original packaging for later use; this can come in handy for transportation purposes.

7.2 Basic equipment

KARL STORZ Set no.: UP 501 S1

- 1 S-PILOT™ UP 501
- 1 Power cord 400 A
- 5 Tubing sets for the suctioning of smoke, gas and fluid, for single use, sterile 031447-03
- 1 SCB connecting cable 20 0901 70, length 100 cm
- 1 One-pedal footswitch 20 0141 30
- 1 Instruction manual

**7.3 Installing and connecting the
device**

i **NOTE:** The S-PILOT™ including the accessories connected may only be used in medical rooms in which the electrical equipment has been installed in accordance with applicable national regulations.

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1. Place the device on a flat surface.


WARNING: The device must be set up in such a way that the power cord can be simply pulled out at any time.

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7.3.1 Connecting the ground line

WARNING: The device is equipped with a connector for attaching a ground line. Carry out potential equalization as per the applicable national regulations.

2. Connect the potential equalization cable to the potential equalization connector.
3. Where necessary, have grounding performed by qualified personnel.

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7.3.2 Connecting the power cord

WARNING: The device may only be operated at the voltage stated on the rating plate.

4. Push the power cord into the power socket as far as it will go.
- ➔ The operating indicator light (status LED) of the On/Off switch lights up orange (standby).

WARNING: Only connect/disconnect the power plug to/from the power supply outside of areas subject to explosion hazards.

WARNING: The instrument may only be operated with the power cord delivered by KARL STORZ or a similar power cord which has a national inspection seal.

WARNING: The line voltage to the device is only definitively disconnected once the power plug has been unplugged.

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**7.3.3 Connecting the footswitch for
direct control**

Connect the single-pedal footswitch 200141 30 to the socket ④, while aligning the red point on the connector of the connecting cable with the red dot on the footswitch socket.

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**7.3.4 Combining the HF device with the
S-PILOT™ in the OR1™ system
via the KARL STORZ-SCB –
automatic smoke evacuation**

i **NOTE:** To prevent the SCB connecting cable from being pulled out accidentally, the SCB connector possesses a protection device.

Pull back the protection device of the SCB connector and insert the connector into one of the SCB sockets. Connect the other end of the cord to the KARL STORZ-SCB (KARL STORZ Communication Bus) control unit or other SCB devices (see KARL STORZ-SCB control NEO System Instruction Manual).

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**7.3.5 Combining the HF device
with the S-PILOT™ via the
KARL STORZ-SCB – automatic
smoke evacuation as
a standalone solution**

• **AUTOCON® II 400**

Directly connect the SCB connectors ⑦ of the S-PILOT™ and AUTOCON® II 400 using the SCB connecting cable 20 0901 70 (see note in section 7.3.4).

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• **ERBE Vio 300D**

In order to connect the ERBE HF device to the S-PILOT™, you require the KARL STORZ-SCB/ERBE Vio interface device 20 0900 13.

1. Connect the SCB connectors ⑦ of the S-PILOT™ and interface device SCB/ERBE Vio 20 0913 32 using the SCB connecting cable 20 0901 70 (see Fig. in section 7.3.5).
2. Connect the interface device SCB/ERBE Vio using the ECB connecting cable 20 0913 70 with the ERBE Vio 300D.

• **KLS Martin maxium ME402**

To connect the KLS Martin maxium HF device to the S-PILOT™ you require the set for SCB/KLS MARTIN maxium 20 0912 90 and, if necessary, a carrier device for SCB/RS-232 interface modules (e.g. KARL STORZ-SCB Interface control 20 0901 20).

1. Connect the SCB connectors ⑦ of the S-PILOT™ and SCB Interface control 20 0901 20 using an SCB connecting cable (see Fig. in section 7.3.5).
2. Connect the SCB Interface control 20 0901 20 using the RS232 connecting cable 20 0918 70 with the KLS Martin maxium.

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7.3.6 Combining the HF device with the S-PILOT™ via the connecting cable UP 004/UP 005 – automatic smoke evacuation as a standalone solution

• **AUTOCON® III 400**

Connect the remote control connection ⑤ of the S-PILOT™ with the corresponding connection of the AUTOCON® III 400 HF device using the connecting cable UP 004 (dia. 3.5 mm).

• **ConMed System 2450 and 5000**

Connect the remote control connection ⑤ of the S-PILOT™ with the corresponding connection of the HF device using the connecting cable UP 004 (dia. 3.5 mm).

• **COVIDIEN Force Triad and Valleylab Force FX**

Connect the remote control connection ⑤ of the S-PILOT™ (dia. 3.5 mm) with the corresponding connection of the HF device (dia. 2.5 mm) using the connecting cable UP 005.

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
7.3.7 Combining the UNIDRIVE® S III with the S-PILOT™ – suction of fluids

Connect the control cable 27 6100 70 between the UNIDRIVE® S III and S-PILOT™ (footswitch socket ④).

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7.3.8 Combining the CALCUSON with the S-PILOT™ – suction of fluids

Connect the control cable 20 7010 70 between the CALCUSON and S-PILOT™ (footswitch socket ④).

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7.4 Commissioning

1. Press the On/Off switch ① (see Fig. A).

The status LED of the On/Off switch lights up green. The tappet of the tube squeeze valve moves backwards and then forwards. Once the test has successfully completed, an acoustic availability signal is emitted (see Section 9.1 – Acoustic signals).

NOTE: If the status LED of the On/Off switch ① is continuously illuminated and the valve does not open and close when the device is switched on, an SD card is probably still in the slot. Please remove it, the system then starts normally.

7.4.1 Inserting the suction tube set

NOTE: The tubing set which has been designated as suitable by KARL STORZ (see Section 10.2) ensures optimal functioning of the S-PILOT™ UP 501.

WARNING: No responsibility for the safe operation of the device can be accepted if any tubing systems other than those specified are used.

WARNING: When inserting the tubing set, it must be ensured that the tube is completely inserted.

WARNING: Before you start the intervention, make sure that the tubing set is not leaking or kinked and that the ends of the tubes are properly secured.

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Opening the tube squeeze valve

i **NOTE:** The tube squeeze valve ③ is closed when not activated.

2. Press the tube squeeze valve button ② (see Fig. B – page 21).

☛ The tube squeeze valve ③ opens (the tappet moves backwards).

i **NOTE:** Every time the button ② is pressed, the tube squeeze valve switches between 'open' and 'closed' ('Toggle' function).

3. Insert the suction tube into the tube squeeze valve (see Fig. C – page 21).

4. Connect one end of the suction tube (marked white) to the suction canister/system (see Fig. E). The connection point on the lid (in the case of the mtp system) is marked red. Ensure that the tube end of the Y-adaptor which is not required is closed (see Fig. D).

i **NOTE:** To prevent the cavity from collapsing, the suction flow can be adjusted at any time via the roller clamp attached close to the patient.

i **NOTE:** Fluid and smoke gas can be suctioned simultaneously via the Y-adaptor.

5. Suction tube 031447-01:
Connect the other end of the suction tube to the trocar (via LUER lock; (see Fig. F)).

Suction tube 031457-01:
Connect the other end of the suction tube to the shaver handpiece (via silicone tube end; (see Fig. G)).

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7.4.2 Test for proper functioning

WARNING: Check the device combination prior to each use to ensure it is functioning correctly. The device/system should not be used if any damage is evident.

WARNING: The underpressure should not be set greater than required for the intervention (max. -0.3 bar). This can be carried out by setting the vacuum pump or setting the pressure regulator in the case of a central vacuum supply.

NOTE: When the device function is interrupted and the squeeze valve cannot be closed, close the shutoff valve on the tube connection of the instrument.

S-PILOT™ with UNIDRIVE® S III

UNIDRIVE® S III: Connect to the power supply. Connect the footswitch 20 0162 30 to the footswitch socket. Connect the motor handpiece. Switch on the device. Wait until the boot sequence is complete.

S-PILOT™: Connect to the power supply.

Connect the control cable 27 6100 70 between the UNIDRIVE® S III and the S-PILOT™.

Following the boot sequence – actuate the first stage of the footswitch.

- ☛ The tube squeeze valve must be open as long as the activation signal is present.

S-PILOT™ with CALCUSON

CALCUSON: Connect to the power supply. Connect the footswitch 20 0142 30 to the footswitch socket. Connect the ultrasonic transducer. Fasten the ultrasonic probe to the connector on the transducer.

S-PILOT™: Connect to the power supply.

Connect the control cable 20 7010 70 between the CALCUSON and S-PILOT™.

Switch on the CALCUSON and S-PILOT™.

- ☛ The pilot lamps ⑤ and ⑦ on the CALCUSON must light up green.

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Actuation of the first stage of the two-stage footswitch 20 0142 30 starts suction (opens the tube squeeze valve), the second stage activates the CALCUSON.

- ☛ The tube squeeze valve must be open as long as the activation signal is present.
- ☛ If it operates properly the green pilot lamp ② lights up when ultrasonic energy is emitted.

S-PILOT™ with AUTOCON® II 400

AUTOCON® II 400: Connect to the power supply. Switch on the device. Wait until the boot sequence is complete.

S-PILOT™: Connect to the power supply.

Connect the SCB connecting cable between the AUTOCON® II 400 and S-PILOT™.

Following the boot sequence – actuate the HF device (cutting/coagulation).

- ☛ The tube squeeze valve must open as long as the activation signal is present (plus shut-off delay of 3 seconds).

S-PILOT™ with AUTOCON® III 400

AUTOCON® III 400: Connect to the power supply. Switch on the device. Wait until the boot sequence is complete.

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Following the boot sequence – actuate the HF device (cutting/coagulation).

- ☛ The tube squeeze valve must open as long as the activation signal is present (plus shut-off delay of 3 seconds).

7.4.3 Decommissioning

- Remove the used suction tube from the device and dispose of it.
- Disconnect the device from the power supply and perform a wipe-down disinfection of the device (see section 'Reprocessing').

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8 Maintenance

8.1 Reprocessing

- UP 501 – S-PILOT™ Pinch valve
- Suction tube set for single use, sterile

General warnings

WARNING: Risk of infection: Incorrectly reprocessed medical devices expose patients, users and third parties to a risk of infection as well as the risk that the medical device may malfunction. Observe the 'Cleaning, Disinfection, Care, and Sterilization of KARL STORZ Instruments' instructions and the accompanying documentation.

WARNING: When carrying out any work on contaminated medical devices, the guidelines of the Employer's Liability Insurance Association and equivalent organizations striving to ensure personal safety must be observed.

WARNING: When preparing and using the solutions, follow the chemical manufacturer's specifications, paying close attention to proper concentration, exposure time and service life. Incorrect concentration may result in damage. Bear in mind the microbiological range of action of the chemicals used.

WARNING: National laws and regulations must be observed.

NOTE: The manual 'Cleaning, Disinfection, Care, and Sterilization of KARL STORZ Instruments' can be downloaded or requested by visiting www.karlstorz.com.

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8.1.1 Reprocessing of the UP 501

WARNING: Always disconnect the device from the mains before cleaning!

WARNING: Liquid must not be allowed to enter the housing under any circumstances.

Manual wipe-down disinfection

The exterior surfaces of the medical device must be wiped clean with a disposable cloth moistened with disinfectant or a ready to use soaked disinfectant cloth. Due to their protein-fixating effect and material incompatibility, alcohol-based agents must not be used. The chemical manufacturer's specifications concerning material compatibility must be observed. At the end of the necessary exposure time for the disinfectant as specified by the manufacturer, wipe the surface with a dry and low lint cloth.

8.1.2 Reprocessing of suction tube set for single use, sterile

WARNING: Suction tube sets for single use must not be reprocessed.

Sterile disposable products are not approved for reprocessing.

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
8.2 Maintenance and safety check

8.2.1 Maintenance

Preventive maintenance is not essential. Regular maintenance can, however, contribute to identifying potential problems before they become serious, thus enhancing the device's reliability and extending its useful service life. Information regarding maintenance services can be obtained from your local representative or from the manufacturer.

8.2.2 Safety checks/repeat inspections as per IEC 62353

WARNING: *Regardless of the accident prevention regulations or testing intervals for medical devices prescribed in different countries, with this device safety checks/repeat inspections as defined by IEC 62353 must be performed and recorded by a qualified electrician at least once a year.*

Visual inspection

1. Check the device and accessories for any mechanical damage which may impair functionality.
2. Check that inscriptions relevant to safety are legible.

Electric measurements

- Measure protective ground resistance in accordance with IEC 62353: The limit values can be found in the current standard.
- Measure earth leakage current in accordance with IEC 62353: The limit values can be found in the current standard.

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Test for proper functioning

Carry out the functional check as per the instruction manual (see section 7.4.2).

Documentation

The safety check/repeat inspection and results must be documented.

WARNING: If any defects or shortcomings which could endanger patients, personnel or others are ascertained during these safety checks/repeat inspections, the device must not be operated until such time as these defects or shortcomings have been eliminated by qualified technical servicing.

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8.3 Servicing and repair

WARNING: Defective items of equipment must be serviced and repaired exclusively by KARL STORZ or persons authorized by KARL STORZ; all repair work must employ original KARL STORZ parts only.

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8.4 Disposal

This device has been marked in accordance with the European Directive on Waste Electrical and Electronic Equipment (WEEE).

WARNING: At the end of its service life, dispose of the device as electronic waste.

Please ask either KARL STORZ GmbH & Co. KG, a KARL STORZ subsidiary or your specialist dealer for information on your local collection point. Within the scope of application of this Directive, KARL STORZ GmbH & Co. KG is responsible for the proper disposal of this device.

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Servicehotline: +49 7461/708 980
E-Mail: technicalsupport@karlstorz.com

8.5 Repair program

During the repair period, you will generally be loaned a device, which must then be returned to KARL STORZ as soon as the repaired device is returned to you.

In Germany you can refer repairs direct to:

KARL STORZ GmbH & Co. KG
Abt. Reparaturservice
Take-off Gewerbepark 83
78579 Neuhausen

Service hotline: +49 7461/708 980
E-mail: technicalsupport@karlstorz.com

In other countries please contact your local KARL STORZ branch or authorized dealer.

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8.6 Important information

In order to prevent infection, it is strictly forbidden to ship contaminated medical devices. All medical devices must be decontaminated on site, to avoid contact and aerogenous infections (among staff). We reserve the right to return contaminated instruments/equipment to the sender.

Repairs, modifications, or expansions which are not carried out by KARL STORZ or by experts authorized by KARL STORZ will invalidate all warranty rights. KARL STORZ gives no guarantee of the correct functioning of devices or instruments which have been repaired by unauthorized third parties.

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8.7 Limitation of liability

KARL STORZ shall be liable for the safe operation, operational reliability, and performance of this equipment only subject to the following conditions:

- all assembly operations, system expansions, readjustments, modifications, or repairs have been performed by a person or persons duly authorized by KARL STORZ;
- the electrical installation of the location in which the device is connected and operated conforms to the applicable statutes and standards, and
- the device has been used in accordance with its operating instructions at all times.

8.8 Warranty

The guarantees provided can be found in the Standard Conditions of Business of KARL STORZ.

The medical device must always be sent to your local subsidiary (see "Subsidiaries" section), even during the warranty period.

Opening the equipment or performance of any repairs or modifications to the equipment by unauthorized persons shall relieve us of any liability for its performance. Any such repair or modification performed during the warranty period shall void all warranty.

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


Availability signal ('Ready')

This signal is issued following a successful self-test.

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Service kontaktieren.

9.2 Troubleshooting

WARNING: Always unplug the device before performing any maintenance work on it!

Error description:

Complete failure of the device (status LED of ON/ Off switch does not light up).

Possible causes:

1. Power cord connector is not properly connected to device socket.
2. Power supply failure.
3. Defective internal device fuse.

Action:

1. Firmly push power supply connector into socket on device.
2. Check that there is electricity to the wall outlet.
3. Contact Service.

WARNING: The user can not replace a defective internal device fuse.

Error description:

Inadequate suction.

Possible causes:

1. Leakage in tubing system.
2. Tubing system blocked.
3. Suction container full.

Action:

1. Check all suction lines, and replace any exhibiting evidence of leakage. Check the seating of the sealing cap.
2. Remove all blockages, and replace any exhibiting evidence of leakage.
3. Empty or replace suction container.

Symptom:

Error during Power on Test (Chap. 6.8.1)

- System error:
Status LED and buzzer: T_{on}/T_{off} 1 s/1 s
- Motor driver:
Status LED and buzzer: T_{on}/T_{off} 1 s/3 s

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Symptom:

SCB connection lost during operation:
Acoustic information signal every 10 s

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- Type of protection against electric shocks: Protection Class I
- Degree of protection against electric shocks: Applied part type CF
- Type of IP protection: IP21

According to IEC 60601-1-2, FCC-B:

Observe the information on electromagnetic compatibility in chapter 11 (pages 38-51).

**9.3.2 Directive compliance
(for UP 501)**

According to the Medical Device Directive (MDD):

This medical product belongs to Class I.

This medical product bears the CE mark in accordance with the Medical Device Directive (MDD) 93/42/EEC.

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9.4 Technical documentation

On request, the manufacturer will provide available circuit diagrams, detailed parts lists, descriptions, instructions for adjustments, and other documentation deemed useful for suitably qualified user personnel authorized by the manufacturer when repairing device components that have been designated as repairable by the manufacturer.

Supply of such technical documentation relating to the device shall not be construed as constituting manufacturer's authorization of user's personnel, regardless of their level of technical training, to open or repair the device.

Explicitly exempted herefrom are those maintenance and repair operations described in this manual.

i NOTE: We reserve the right to make engineering modifications in the interest of promoting technological progress and generating performance improvements without obligation on our part to submit prior notice thereof.

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**10 Spare parts,
recommended accessories**
10.1 Spare parts/accessories

Article	Cat. no.
One-pedal footswitch, one-stage	20 0141 30
Control cable (connecting cable between the UNIDRIVE® S III and the S-PILOT™), length 42 cm	27 6100 70
SCB connecting cable, length 100 cm	20 0901 70
Power cord (grounded)	400 A
Power cord 'Hospital Grade' (USA)	400 B
Instruction manual	96116047 D

10.2 Accessories

Article	Cat no.
Control cable (connecting cable between the UNIDRIVE® S III / CALCUSON and the S-PILOT™), length 100 cm	20 7010 70
Single-use tubing set for the suctioning of smoke, gas and fluid. Sterile, 10 per pack.	031447-10
Single-use tubing set for suction, sterile, 10 per pack.	031457-10
VACUsafe suction bag 2 l with disposable filter. Color code: green. For use with VACUsafe canister 2 l. 48 per pack.	030220-48
VACUsafe canister 2 l. For use with the two VACUsafe suction bags 2 l. 18 per pack.	030020-18
VACUsafe connecting tube, 30 cm, with green multiadaptors. Non-sterile, 10 per pack.	030648-10
Smoke gas filter for single use. Non-sterile, 10 per pack.	031111-10
S-PILOT™ connecting cable, dia. 3.5 mm, length 300 cm, for use with AUTOCON® III 400, ConMed System 2450 or 5000	UP 004
S-PILOT™ connecting cable, dia. 2.5 mm, length 300 cm, for use with Covidien Force Triad or Covidien Valleylab Force FX	UP 005
Set for SCB/ERBE VIO Interface Box	20 0900 13
Set for SCB/KLS MARTIN maxium	20 0912 90
KARL STORZ SCB Interface control Set	20 0905 01

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**11 Electromagnetic
Compatibility (EMC)
information**

WARNING: Medical electrical devices are subject to special precautions regarding electromagnetic compatibility (EMC). Observe the EMC instructions in this section during installation and operation.

The KARL STORZ S-PILOT™ model UP 501 corresponds to IEC 60601-1-2 [CISPR 11 Class B] and therefore meets the EMC requirements of the Medical Device Directive (MDD) 93/42/EEC.

The limits used are designed to provide reasonable protection against the typical electromagnetic interference to be expected in a medical environment. The KARL STORZ S-PILOT™ model UP 501 is a Group 1 device (as per CISPR 11).

Group 1 contains all of the 'equipment and systems which generate or use RF energy only for their internal functioning'.

NOTE: The tables and guidelines that are included in this section provide information to the customer or user that is essential in determining the suitability of the device or system for the electromagnetic environment of use, and in managing the electromagnetic environment of use to permit the device or system to perform its intended use without causing interferences in other devices or systems or non-medical electrical devices. If this device does cause electromagnetic interferences with other devices, the user is encouraged to try to correct the interferences by one or more of the following measures:

- reorient or relocate the receiving device
- increase the separation between the individual devices
- connect the devices to different electrical circuits.

If you have any further questions, please contact your local representative or our service department.

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Electromagnetic Compatibility (EMC) information

 **WARNING:** Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.

 **NOTE:** The device has been tested in connection with the following devices and can be operated in their immediate vicinity or stacked with these.
ENDOFLATOR® 40 – UI 400
ENDOFLATOR® 50 – UI 500
Electronic ENDOFLATOR® – 26 4305 20-1
THERMOFLATOR® – 26 4320 20-1
UNIDRIVE® S III – 20 7010 20-1

 **WARNING:** Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the S-PILOT™ UP 501, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

 **WARNING:** Use of accessories, transducers and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.

 **NOTE:** The device S-PILOT™ tube squeeze valve does not feature any ESSENTIAL PERFORMANCE CHARACTERISTICS with regard to its function and has been exclusively tested regarding immunity to electromagnetic interferences to the effect that the device does not fail and internal assemblies are not damaged.

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Cables which have been shown to comply with IEC 60601-1-2:				
Cable type	Shielded	Length [m]	Ferrite	Used for
PA	No	> 3	No	Potential equalization
SCB Cable	Yes	≥ 1	No	Connecting SCB network, other SCB units
20 0141 30	Yes	3	No	One-pedal-Footswitch
27 6100 70	Yes	0.42	No	Control cable (connecting cable between the S-PILOT™ and the UNIDRIVE® S III)
20 7010 70	Yes	1	No	Control cable (connecting cable between the S-PILOT™ and UNIDRIVE® S III /CALCUSON)
UP 004	No	3	No	Connecting cable between the S-PILOT™ and AUTOCON® III 400/ConMed system 2450 or 5000
UP 005	No	3	No	Connecting cable between the S-PILOT™ and Covidien Force Triad or Covidien Valleylab Force FX
Power Cord	No	3	No	Power supply

Table 1 Guidance and manufacturer's declaration – electromagnetic emissions		
The S-PILOT™ model UP 501 is intended for use in the electromagnetic environment specified below. The customer or user of the S-PILOT™ model UP 501 should ensure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment – guidance
RF emissions CISPR 11	Group 1	The S-PILOT™ model UP 501 uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The S-PILOT™ model UP 501 is suitable for use in all establishments including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/flicker emissions IEC 61000-3-3	Complies	

Table 2			
Guidance and manufacturer's declaration – electromagnetic immunity			
The S-PILOT™ model UP 501 is intended for use in the electromagnetic environment specified below. The customer or user of the S-PILOT™ model UP 501 should ensure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic Discharge (ESD) IEC 61000-4-2	± 6 kV contact ± 8 kV air	± 6 kV contact ± 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines	± 2 kV for power supply lines ± 1 kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5 % U_T^* (>95 % dip in U_T) for 0.5 cycle 40 % U_T (60 % dip in U_T) for 5 cycles 70 % U_T (30 % dip in U_T) for 25 cycles <5 % U_T (>95 % dip in U_T) for 5 sec	<5 % U_T^* (>95 % dip in U_T) for 0.5 cycle 40 % U_T (60 % dip in U_T) for 5 cycles 70 % U_T (30 % dip in U_T) for 25 cycles <5 % U_T (>95 % dip in U_T) for 5 sec	Mains power quality should be that of a typical commercial or hospital environment. If the user of the equipment or system requires continued operation during power mains interruptions, it is recommended that the equipment or system be powered from an uninterruptible power supply.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
* Note: U_T is the a.c. mains voltage prior to application of the test level.			

<p align="center">Table 4 Guidance and manufacturer's declaration – electromagnetic immunity – for equipment and systems that are not life-supporting</p>			
<p align="center">The S-PILOT™ model UP 501 is intended for use in the electromagnetic environment specified below. The customer or user of the S-PILOT™ model UP 501 should ensure that it is used in such an environment.</p>			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
<p>Conducted RF IEC 61000-4-6</p>	<p>3 V_{rms} 150 kHz to 80 MHz</p>	<p>3 V_{rms}</p>	<p>Portable and mobile RF communications equipment should be used no closer to any part of the S-PILOT™ model UP 501, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance: $d = [3.5/3]\sqrt{P}$ $d = [3.5/3]\sqrt{P}$ 80 MHz to 800 MHz $d = [7/3]\sqrt{P}$ 800 MHz to 2.5 GHz</p> <p>Where P is the maximum output power rating of the transmitter in watts [W] according to the transmitter manufacturer and d is the recommended separation in meters [m].</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey^a, should be less than the compliance level in each frequency range ^b.</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol: </p>
<p>Radiated RF IEC 61000-4-3</p>	<p>3 V/m 80 MHz to 2.5 GHz</p>	<p>3 V/m</p>	
<p>Note 1: At 80 MHz and 800 MHz, the higher frequency range applies. Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection by structures, objects and people.</p>			
<p>a Field strengths from fixed transmitters, such as base stations for radio [cellular/cordless] telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the unit is used exceeds the applicable RF compliance level above, the unit should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the unit.</p> <p>b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.</p>			

Table 6

Recommended separation distances between portable and mobile RF communications equipment and the S-PILOT™ model UP 501

The S-PILOT™ model UP 501 is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or user of the unit can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the unit as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter [W]	Separation distance d [m] according to frequency of transmitter		
	150 kHz to 80 MHz $d = \left[\frac{3.5}{3} \right] \sqrt{P}$	80 MHz to 800 MHz $d = \left[\frac{3.5}{3} \right] \sqrt{P}$	800 MHz to 2.5 GHz $d = \left[\frac{7}{3} \right] \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance **d** in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where **P** is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

Note 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

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KARL STORZ Endoscopy Canada Ltd.
7171 Millcreek Drive, Mississauga
Ontario, L5N 3R3, Canada
Phone: +1 905 816-4500, Fax: +1 905 858-4599
Toll free phone: 1-800-268-4880 (Canada only)
Toll free fax: 1-800-482-4198 (Canada only)
E-Mail: info.canada@karlstorz.com

KARL STORZ Endoscopy-America, Inc.
2151 East Grand Avenue
El Segundo, CA 90245-5017, USA
Phone: +1 424 218-8100, Fax: +1 424 218-8526
Toll free phone: 800 421-0837 (USA only)
Toll free fax: 800 321-1304 (USA only)
E-Mail: info@kseam.com

KARL STORZ Veterinary Endoscopy-America, Inc.
175 Cremona Drive
Goleta, CA 93117, USA
Phone: +1 805 968-7776, Fax: +1 805 685-2588
E-Mail: info@karlstorvet.com

KARL STORZ Endoscopy Latino-America, Inc.
815 N. W. 57th Avenue, Suite 480
Miami, FL 33126-2042, USA
Phone: +1 305 262-8980, Fax: +1 305 262-8986
E-Mail: info@ksela.com

KARL STORZ Endoscopia
Miramar Trade Center,
Edificio Jerusalem, Oficina 108,
La Habana, Cuba
Phone: +53 72041097, Fax: +53 72041098
E-Mail: karlstorz@enet.cu

KARL STORZ Endoscopia México S.A de C.V.
Lago Constanza No 326,
Col. Granada Del. Miguel Hidalgo
C.P. 11520 México D.F.
Phone: +52 (55) 1101 1520
E-Mail: mx-info@karlstorz.com

KARL STORZ Marketing América Do Sul Ltda.
Rua Joaquim Floriano, nº. 413, 20º andar - Itaim Bibi,
CEP-04534-011 São Paulo, Brasil
Phone: +55 11 3526-4600, Fax: +55 11 3526-4680
E-Mail: info@karlstorz.com.br

KARL STORZ Endoscopia Argentina S.A.
Zufriategui 627 6º Piso, B1638 CAA - Vicente Lopez
Provincia de Buenos Aires, Argentina
Phone: +54 11 4718 0919, Fax: +54 11 4718 2773
E-Mail: info@karlstorz.com.ar

KARL STORZ Endoskopi Norge AS
Østensjøveien 15B
0661 Oslo, Norway
Phone: +47 6380 5600, Fax: +47 6380 5601
post@karlstorz.no

KARL STORZ Endoskop Sverige AB
Storsåtragränd 14
12739 Skärholmen, Sweden
14108 Kungens Kurva, Sweden
Phone: +46 8 505 648 00
E-Mail: kundservice@karlstorz.se

KARL STORZ Endoscopy Suomi OY
Valimotie 1 B, 4. kerros
00380 Helsinki, Finland
Phone: +358 (0)96824774, Fax: +358 (0)968247755
E-Mail: asiakaspalvelu@karlstorz.fi

KARL STORZ GmbH & Co. KG
Representation Office
Kestucio g. 59
08124 Vilnius, Lithuania
Phone: +370 68567000
E-Mail: ausra.kukeniene@karlstorz.com

KARL STORZ Endoskopi Danmark A/S
Skovlytoften 33
2840 Holte, Denmark
Phone: +45 45162600, Fax: +45 45162609
E-Mail: marketing@karlstorz.dk

KARL STORZ Endoscopy (UK) Ltd.
415 Perth Avenue, Slough
Berkshire, SL1 4TQ, United Kingdom
Phone: +44 1753 503500, Fax: +44 1753 578124
E-Mail: info-uk@karlstorz.com

KARL STORZ Endoscopie Nederland B. V.
Displayweg 2
3821 BT Amersfoort, Netherlands
Phone: +31 (0)33 4545890
E-Mail: info-nl@karlstorz.com

KARL STORZ Endoscopy Belgium N. V.
Phone: +32 473 810 451
E-Mail: info@stopleer.be

KARL STORZ Endoscopie France S. A. S.
12, rue Georges Guynemer, Quartier de l'Europe
78280 Guyancourt, France
Phone: +33 1 30484200, Fax: +33 1 30484201
E-Mail: marketing@karlstorz.fr

KARL STORZ Endoskop Austria GmbH
Landstraße Hauptstr. 148/1/G1
1030 Wien, Austria
Phone: +43 1 71 56 0470, Fax: +43 1 71 56 0479
E-Mail: storz-austria@karlstorz.at

KARL STORZ Endoscopia Ibérica S. A.
Parque Empresarial San Fernando
Edificio Munich - Planta Baja
28830 Madrid, Spain
Phone: +34 91 6771051, Fax: +34 91 6772981
E-Mail: info-es@karlstorz.com

KARL STORZ Endoscopy Italia S. r. l.
Via dell'Artigianato, 3
37135 Verona, Italy
Phone: +39 045 8222000, Fax: +39 045 8222001
E-Mail: info-ita@karlstorz.com

KARL STORZ Adria Eos d.o.o.
Capraška 6
10000 Zagreb, Croatia
Phone: +385 1 6406 070, Fax: +385 1 6406 077
E-Mail: info@karlstorz.hr

KARL STORZ Endoskopija d.o.o.
Verovškova c. 60A
1000 Ljubljana, Slovenia
Phone: +385 1 620 5880, Fax: +385 1 620 5882
E-Mail: pisarna@karlstorz.si

KARL STORZ Endoskop Polska Marketing Sp. z o.o.
ul. Bojkowska 47
44-100 Gliwice, Poland
Phone: +48 32 706 13 00, Fax: +48 32 706 13 07
E-Mail: info-pl@karlstorz.com

KARL STORZ Endoszkóp Magyarorszá g Kft.
Toberek utca 2. fsz. 17/b
HU-1112 Budapest, Hungary
Phone: +36 195 096 31, Fax: +36 195 096 31
E-Mail: info-hu@karlstorz.com

KARL STORZ Endoscopia Romania srl
Str. Prof. Dr. Anton Colorian, nr. 74, Sector 4
041393 Bukarest, Romania
Phone: +40 (0)31 4250800, Fax: +40 (0)31 4250801
E-Mail: info-ro@karlstorz.com

KARL STORZ Endoskope Greece E.P.E.*
Patriarhou Grigoriou E' 34
54248 Thessaloniki, Greece
Phone: +30 2310 304868, Fax: +30 2310 304862
E-Mail: info@karlstorz.gr

*Repair & Service Subsidiary

KARL STORZ Industrial**
Gedik Is Merkezi B Blok
Kat 5, D 38-39, Bagdat Cad. No: 162
Maltepe Istanbul, Turkey
Phone: +90 216 442 9500, Fax: +90 216 442 9030

**Sales for Industrial Endoscopy

000 KARL STORZ Endoscopy - WOSTOK
Derbenyevskaya nab. 7, building 4
115114 Moscow, Russia
Phone: +7 495 983 02 40, Fax: +7 495 983 02 41
E-Mail: kste-wostok@karlstorz.ru

TOV KARL STORZ Ukraine
Obolonska naberezhna, 15
building 3, office 3
04210 Kiev, Ukraine
Phone: +380 44 42668-14, -15, -19
Fax: +380 44 42668-20
E-Mail: marketing@karlstorz.com.ua

KARL STORZ GmbH & Co. KG
Representation Office
Sabit Orudschow 1184, apt. 23
1025 Baku, Azerbaijan
Phone: +99 450 613 30 60
E-Mail: taley.bashirov@karlstorz.com

KARL STORZ ENDOSKOPE
East Mediterranean and Gulf S.A.L.
Block M, 3rd Floor
Beirut Souks, Weyqand Street
2012 3301 Beirut, Lebanon
Phone: +961 1 999390, Fax: +961 1 999391
E-Mail: info@karlstorz-emg.com

KARL STORZ Endoscopy (South Africa) (Pty) Ltd.
P.O. 6061, Roggebaai 8012
Cape Town, South Africa
Phone: +27 21 417 2600, Fax: +27 21 421 5103
E-Mail: info@karlstorz.co.za

TOO KARL STORZ Endoskopij Kasachstan
6. Saryarka str. BC „Arman“, office314
010000 Astana, Republic of Kazakhstan
Phone: +7 7172 552-549, 552-788, Fax: -444
E-Mail: info@karlstorz.kz

KARL STORZ ENDOSKOPE
East Mediterranean & Gulf (branch)
Building West Side 7A - Unit 7WA - 3008
Dubai Airport Free Zone, P.O. Box 54983
Dubai - United Arab Emirates
Phone: +971 (0)4 2958887, Fax: +971 (0)4 3205282
Service Hotline: +971 (0)4 3415882
E-Mail: info-gne@karlstorz-emg.com

KARL STORZ Endoscopy India Private Ltd.
D-181, Okhla Industrial Area, Phase-1
New Delhi 110020, India
Phone: +91 11 43743000-30,
Fax: +91 11 268129 86
E-Mail: karlstorz@vsnl.com

KARL STORZ GmbH & Co. KG
Resident Representative Office
14th Floor, Kumho Asiana Plaza
39 Le Duan, District 1
Ho Chi Minh City, Vietnam
Phone: +84 8 3823 8000, Fax: +84 8 3823 8039
E-Mail: infovietnam@karlstorz.com

KARL STORZ Endoscopy Australia Pty. Ltd.
15 Orion Road Lane Cove NSW 2066
P O Box 50 Lane Cove NSW 1595, Australia
Phone: +61 (0)2 9490 6700
Toll free: 1800 996 562 (Australia only)
Fax: +61 (0)2 9420 0695
E-Mail: info@karlstorz.au

KARL STORZ Endoscopy China Ltd., Hong Kong
Room 2503-05, 25F AXA Tower, Landmark East,
No. 100 How Ming Street, Kwun Tong, Kowloon,
Hong Kong, People's Republic of China
Phone: +852 28 65 2411, Fax: +852 28 65 4114
E-Mail: inquiry@karlstorz.com.hk

KARL STORZ Endoscopy (Shanghai) Ltd.
Beijing Branch,
Room 1805-1807, Building B, 18F Beijing IFC
No. 8, Jianguomenwai Street, Chaoyang District,
Beijing, 100022, People's Republic of China
Phone: +86 10 5638188, Fax: +86 10 5638199
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Shanghai Branch
36F, Tower 1 Grand Gateway
No. 1 Hongquiao Road, Xuhui District,
Shanghai, 200030, People's Republic of China
Phone: +86 21 6113-1188, Fax: +86 21 6113-1199
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Chengdu Branch
Room 805, 8F Jin Jiang International Building
No. 1 West Linjiang Road, Wuhou District,
Chengdu, 6100414, People's Republic of China
Phone: +86 28 86587977, Fax: +86 28 86587975
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Shenyang Branch,
Room 2001-2005, 20F N-MEDIA International Center,
No. 167 Youth Avenue, Shenhe District,
Shenyang, 110014, People's Republic of China
Phone: +86 24 23181118, Fax: +86 24 23181119
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Guangzhou Branch
Room 02B, 35F Teem Tower,
No. 208 Tianhe Road, Tianhe District,
Guangzhou, 510620, People's Republic of China
Phone: +86 20 8732-1281, Fax: +86 20 8732-1286
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy Asia Marketing Pte Ltd.
No. 8 Commonwealth Lane #03-02
Singapore 149555, Singapore
Phone: +65 69229150, Fax: +65 69229155
E-Mail: infoasia@karlstorz.com.sg

KARL STORZ Endoscopy Singapore Sales Pte Ltd
No. 8 Commonwealth Lane #03-02
Singapore 149555, Singapore
Phone: +65 69229150, Fax: +65 69229155
E-Mail: infoasia@karlstorz.com.sg

KARL STORZ Endoscopy Korea Co. Ltd.
I-Park Bundang1, 102dong, 8F
239 Jeongjail-Ro, Bundang-Gu
Seongnam, 463-859, Korea
Phone: +82-70-4350-7474, Fax: +82-70-8277-3299
E-Mail: stephen.Lee@karlstorz.com

KARL STORZ Endoscopy Taiwan Ltd.
12F, No. 192, Sec. 2, Chung Hsin Rd.,
Sindian District, New Taipei City, Taiwan
Phone: +886 933 014 160, Fax: +886 2 8672 6399
E-Mail: mingwang@karlstorz.tw

KARL STORZ Endoscopy Japan K. K.
Stage Bldg. 8F, 2-7-2 Fujimi
Chiyoda-ku, Tokyo 102-0071, Japan
Phone: +81 3 6380-8622, Fax: +81 3 6380-8633
E-Mail: info@karlstorz.co.jp

www.karlstorz.com

STORZ
KARL STORZ—ENDOSKOPE

KARL STORZ GmbH & Co. KG

Mittelstraße 8
78532 Tuttlingen

Postfach 230
78503 Tuttlingen
Germany

Telefon: +49 (0)7461 708-0
Telefax: +49 (0)7461 708-105
E-Mail: info@karlstorz.com
Web: www.karlstorz.com