STORZ

KARL STORZ — ENDOSKOPE

GEBRAUCHSANWEISUNG
TP101 – TELE PACK +

INSTRUCTION MANUAL
TP101 – TELE PACK +

MANUAL DE INSTRUCCIONES
TP101 – TELE PACK +

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Thank you for your expression of confidence in the KARL STORZ brand name. Like all of our other products, this product is the result of years of experience and great care in manufacture. You and your organization have decided in favor of a modern, high-quality piece of equipment from KARL STORZ.

This manual is intended to serve as an aid in the proper installation, connection, and operation of the TELE PACK +. All the necessary details and actions are clearly explained. We therefore ask that you read this manual carefully before proceeding to work with the instrument. Keep this manual available for easy reference.

WARNING: KARL STORZ medical devices are not cleaned and disinfected when delivered. These hygiene measures must be performed before the first use/set-up. Cleaning/disinfection of the surfaces is part of the ongoing disinfection measures.

KARL STORZ is continuously working on the further development of all products. Please appreciate that changes to the scope of supply in form, equipment and technology are possible for this reason. Therefore, no claims may be derived from the information, figures and descriptions in this manual.

- © All product illustrations, product descriptions and texts are the intellectual property of KARL STORZ SE & Co. KG. Their use and reproduction by third parties require the express approval of KARL STORZ SE & Co. KG. All rights reserved.

Rx only Federal (USA) law restricts this device to sale by or on the order of a physician.

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
It is recommended that the suitability of the products for the intended procedure be checked prior to use.

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**2 Controls, displays,
connectors and their uses**

①	Monitor and active touch surface
②	Light port
③	On/Off button
④	Side connection area
⑤	USB 2.0 (2x)
⑥	Connection for stroboscopy microphone
⑦	Camera connection 1 C-Line
⑧	Camera connection 2 X-Line
⑨	Mount for camera head holder
⑩	Carrying handle
⑪	VESA 100 mount
⑫	Cable winding
⑬	Speaker
⑭	Potential equalization
⑮	Power cord socket
⑯	Speaker/headphone connection
⑰	Microphone connection
⑱	DVI Out (external monitor)
⑲	USB 2.0 (4x)
⑳	Network socket
㉑	Kensington Lock connection

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3 Symbols employed

	Follow instructions for use
	Standby/On/Off
	Type CF/Defib. resistant applied part according to IEC 60601-1
	Type BF applied part according to IEC 60601-1
	Lateral stroboscopy microphone input
	USB 2.0
	Line out
	Rear microphone input
	Network connection
	Potential equalization connector
	Control of pollution caused by electronic information products (China RoHS)
	This device has been marked in accordance with the European Directive on waste electrical and electronic equipment (WEEE).
	Caution
	Hot surface
	Manufacturer

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	Catalogue number
	Serial number
	Number of products in the product packaging
	Keep dry
	Fragile, handle with care
	Consult instructions for use

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4 Description of the device

The TELE PACK + is a device which combines the following functions in a compact system:

- Illumination
- Image display
- Image processing and
- Comprehensive treatment documentation

Thanks to its novel design, which includes a complete glass front for optimal cleaning and as a touchscreen, it combines state-of-the-art technology with the unique image quality of KARL STORZ camera heads and videoscopes.

The TELE PACK + is compatible with the latest KARL STORZ videoscopes and 1 chip full HD camera heads.

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5 Safety instructions

**5.1 Explanation of Warnings
and Cautions**

The words **Warning**, **Caution** and **Note** convey special meanings. Wherever they are used in this manual, they should be carefully reviewed to ensure the safe and effective operation of this system. To make these words stand out more clearly, they are accompanied by a pictogram.

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WARNING: A warning indicates that the personal safety of the patient or user is at risk. Failure to observe a warning could result in injury to the patient or user.

CAUTION: A Caution indicates that particular service procedures or precautions must be followed to avoid possible damage to the system.

NOTE: A note indicates special information about operating the system, or clarifies important issues.

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WARNING: Read this instruction manual thoroughly before using the equipment. Read the chapter on safety instructions carefully to avoid putting your patients, personnel or yourself at risk.

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5.2 Safety regulations

WARNING: The electrical installations in the operating room in which the unit is connected and operated must comply with the applicable IEC standards.

WARNING: The instruction manuals and interface specifications for medical devices and/or system components used in combination must be observed precisely.

WARNING: Combinations of medical devices are only assured to be safe if

- these are identified as such in the respective instruction manuals or
- the intended use and interface specifications of the devices used in combination permit this.

WARNING: Do not open the device! Danger of an electric shock. Have servicing performed exclusively by the manufacturer or by personnel authorized by the manufacturer. Any opening of the device components by unauthorized personnel will void the warranty.

WARNING: Risk of electric shock! This device must only be connected to a power supply network with a protective conductor.

WARNING: Test this device prior to each surgical procedure to ensure that it functions correctly. The device should not be used if any damage is evident.

WARNING: Keep the device out of reach of patients.

WARNING: For safety reasons, do not simultaneously touch the device output sockets and the patient.

WARNING: Modifications to the device are not permitted.

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-  **WARNING:** The device must only be operated under the ambient conditions specified.
-  **WARNING:** The device may be subject to random failures. Always keep a spare device at hand if the device is essential for your intervention.
-  **WARNING:** Additional devices that are connected to electrical medical devices must conform to the respective IEC or ISO standards (e.g., IEC 60950 for data processing equipment). Furthermore, all configurations must satisfy the requirements for medical electrical systems. Any person who connects additional devices to electrical medical devices configures a medical system, and is therefore responsible for ensuring that this system complies with the requirements for medical electrical systems. It must be taken into account that local regulations take priority over the aforementioned requirements. If in doubt, please contact your local representative or technical support.
-  **WARNING:** Do NOT use the device in the presence of flammable anesthetics.
-  **WARNING:** The device must not be used in oxygenated environments.
-  **WARNING:** Connect or disconnect the power plug to or from the power supply only outside of explosion hazard areas.
-  **WARNING:** The device may only be operated with the power cable delivered by KARL STORZ or a similar power cable which has a national inspection seal.
-  **WARNING:** Make sure to observe national regulations/laws when disposing of the accessories.

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-  **WARNING:** The correct configuration of the device must be checked after each software update.
-  **WARNING:** The device must be set up so as to ensure that the power cord can simply be pulled out at any time.
-  **WARNING:** Grounding reliability can only be achieved when the device is connected to a 'Hospital Only' or 'Hospital Grade' outlet (i.e., approved for use in an operating room environment).
-  **WARNING:** Medical electrical devices are subject to special precautions regarding Electromagnetic Compatibility (EMC). Observe the EMC instructions in chapter 11 during installation and commissioning.
-  **WARNING:** Maintenance may only be performed by specialists authorized by KARL STORZ.
-  **WARNING:** Always unplug the device before carrying out any cleaning and maintenance work.
-  **WARNING:** At the end of its useful operating life, dispose of the device as electronic scrap.
-  **WARNING:** The device is equipped with a connector for potential equalization. Carry out the potential equalization in accordance with the applicable national regulations.
-  **WARNING:** The light source is transmitting high-energy light leading to high surface temperatures on the light cable and on the tip of the endoscope. Tissue damage is possible.
-  **WARNING:** Always use the lowest possible light setting for your intervention.
-  **WARNING:** Danger of glare! Never look into the light outlet when the light system is switched on.
-  **WARNING:** When changing the light cable connector (adaptor) make sure that the light source is switched off. Before replacing the light adaptor, allow the light source to cool down first.

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 **WARNING:** The device screen is made of glass. This glass can break if the device is dropped on a hard surface or receives a substantial impact. If the glass breaks, do not touch the glass parts of the device or attempt to remove the broken glass from the device. Stop using the device until the glass is replaced by qualified service personnel.

 **WARNING:** The device may only be operated within specified ambient conditions.

 **WARNING:** When using high-frequency surgical equipment, keep the working part of the active electrode in the field of view of the operator to avoid accidental high frequency burns. Avoid contact between the active electrode and metal parts of the endoscope.

 **WARNING:** In the case of loss of surgical image for greater than 4 seconds, a power cycle may be required to restore the image.

 **WARNING:** When using two camera heads/videoscopes, always make sure that you have chosen the right source for the current procedure.

 **WARNING:** This medical device provides and manages sensitive information. Therefore, access to the device must be restricted to persons who are authorized to handle patient data.

 **CAUTION:** Regular back-ups must be performed to safeguard data. KARL STORZ is not liable for loss of data. Responsibility for the data lies with the operator.

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- i NOTE:** All data stored internally on the device is stored in encrypted form. The data must be exported onto a USB storage medium before servicing.
- i NOTE:** Upon saving the data to USB, patient data which is legible as plain text is stored on the storage medium.

5.3 Intended use

5.3.1 Intended purpose

The TELE PACK + is a combination device comprising a light source for illumination, camera control unit (CCU) for image processing and documentation as well as a monitor for visualization. It is intended for endoscopic diagnostic and surgical procedures as well as for stroboscopy.

5.3.2 Contraindications

No contraindications relating directly to the medical device are currently known. The responsible physician must decide whether the foreseen application is admissible based on the general condition of the patient.

5.3.3 Patient profile

The use of the TELE PACK + is not limited to a certain patient profile (sex, age, weight etc.). The attending physician must assess a patient's health and declare it adequate.

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5.3.4 User profile (specialist physician)

- Recognized medical skills for endoscopic interventions in the relevant specialization
- Adequate powers of comprehension to assess the situation rationally at all times during the application
- Adequate language skills in at least one of the languages used on the device and in the instruction manual
- Be thoroughly trained in the operation and use of the device
- Knowledge of the contents of the instruction manual
- No physical impairments that could diminish perception of activation and visual message

5.3.5 User profile (assistant)

- Completed training as medical assistant
- Adequate language skills in at least one of the languages used on the device and in the instruction manual
- Be thoroughly trained in the operation and use of the device
- Knowledge of the contents of the instruction manual
- No physical impairments that could diminish perception of activation and visual message

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5.4 Safety precautions at the site of installation

The device may only be used in medical rooms whose electrical systems have been installed in accordance with applicable national regulations. It is not intended for use in potentially explosive atmospheres. This means, among other things:

- Wherever easily combustible and explosive inhalation anesthetics or mixtures thereof are used, the device must not be operated in the demarcated hazard zone.
- This also applies for easily combustible and explosive chemicals, e.g., skin disinfectants and fast-acting surface disinfectants.

The device is equipped with a connector for potential equalization. It should be connected in accordance with the applicable national regulations.

WARNING: The device should be placed on a flat surface, in a video cart or fixed in place. Do not block air inlets/outlets as overheating may result in fire hazard or device failure during procedure.

WARNING: Ensure sufficient air circulation for the device.

WARNING: The device is not intended for use in potentially explosive atmospheres. Do not operate the device within the demarcated hazard zone while explosive anesthetic gases are in use.

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5.5 Safety precautions when operating the device

It is the user's responsibility to make sure the device is safe and operating properly before use.

NOTE: The documentation accompanying the accessories used must also be observed.

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Carefully remove the TELE PACK + and accessories from their packaging. Check the delivery for missing items and evidence of shipping damage.

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Please file any complaints with the manufacturer or supplier immediately. If possible, retain the original packaging materials for later use; these can come in handy if the device has to be transported.

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- 1 TELE PACK + TP101
- 1 Power cord
- 1 Instruction manual

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6.3 Preparing the device for use

6.3.1 Installation

1. Place the device on a level surface. Do not cover the lateral air vents.

NOTE: Alternatively, the device can be mounted on a VESA 100 holding device.

! WARNING: The line voltage to the device is only definitively disconnected once the power plug has been unplugged.

6.3.2 Connecting the potential equalization

NOTE: The device's ground line should be installed by a qualified electrician.

The TELE PACK+ is equipped with a connector for potential equalization.

1. Connect the potential equalization line with the connector for potential equalization on the device.
2. Connect the potential equalization line to the connector of the treatment room.

6.3.3 Connecting USB devices

The USB interfaces of the TELE PACK+ are suitable for connecting up the following peripherals:

- Keyboard/mouse
- Printer
- Footswitch
- USB stick, FAT32 formatted
- USB to ACC adaptor

6.3.4 Connecting the power supply

! WARNING: Only operate the device with the voltage stated on the manufacturer's identification plate.

1. Before connecting the devices to the electrical supply line, check that the values on the power sockets concur with the values on the device's identification plate.
2. Insert the device plug of the power cord into the power supply cord receptacle of the TELE PACK+.
3. Insert the TELE PACK+ power plug into the power sockets.

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6.3.5 Switching on

To switch on the device, press the ON/OFF key on the front of the device and hold down for 2.4 seconds.

6.3.6 Connecting the camera head/ light cable/videoscopes

Insert the camera connector into the corresponding socket. Insert the light cable into the designated connection.

i NOTE: Please observe the compatibility list (chapter 10).

i NOTE: Please note that 2 cameras can be connected to the device at the same time; however, only one camera can be displayed. If you have connected 2 cameras, you can select the relevant image with the 'Swap Cameras' button in the menu (see chapter 7.2.2).

i NOTE: Please note that the two camera connections offer different degrees of protection against electric shocks. The symbols attached are explained in chapter 3.

i NOTE: The use of original KARL STORZ light cables only is recommended. Light cables from other manufacturers may not be optimally adjusted for light transmission.

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6.3.7 White balance

i NOTE: Perform a white balance before each use in order to adjust the perception of white by the camera and the ability to display colors properly.

1. Ensure that the imaging instrument* is connected to the TELE PACK + and is displaying a live image on the monitor.
2. Connect the light source to the imaging instrument via the light cable.
3. Point the tip of the imaging instrument at any clean white surface (e.g., white gauze).
4. Make sure that the light source emits enough light to sufficiently illuminate the white surface. The monitor should display a completely white and well illuminated live image.
5. Press the  button in the live menu to perform the white balance.
6. Look at the monitor to determine whether the white balance was successful or not.

*A camera head with a connected endoscope, videoscope or exoscope

6.3.8 Connecting a microphone

Insert the microphone into the microphone input socket.

i NOTE: A stroboscopy and a normal microphone can be connected simultaneously.

6.3.9 Connecting an external monitor

There is a DVI-D video output on the rear of the device for connecting up an external monitor.

i NOTE: The connected monitor must support a resolution of 1920 x 1080 pixels.

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6.3.10 Attaching the camera head holder TP002

Using the TP002 camera head holder, the camera head can be stored between two interventions, and a light cable can also be attached to the intended receptacle. In the packaging of the camera head holder there is a socket for the camera head holder, the holder itself and 2 screws.

i NOTE: Make sure that the light source is turned off when attaching the light cable to the mount on the camera head holder.

1. Use a screwdriver to remove the screws on the rear left side of the device (see 9).
2. Place the camera head holder mount flush with the device.
3. Use a screwdriver to screw the screws into the two holes to fix the mount.
4. Place the camera head holder on the mount. The holder can be placed in 2 positions to ensure proper function and to stow the holder in case of mobile use of the device for transport.

6.3.11 Function test

1. Check the TELE PACK+ and the camera head or videoscope for external damage.
2. Ensure that the camera head cable and/or the videoscope cable is not damaged or kinked.
3. Start the light source by pressing the  button in the menu.
4. Point the camera head and/or the videoscope towards an object and check the quality of the display.

i NOTE: The device and accessories may only be used once the test for proper functioning has been successfully completed.

6.3.12 Switching off

To switch off the device, hold the ON/OFF key down for 2.4 seconds.

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7 Operating instructions

7.1 Basic information on operation

Operation of the TELE PACK + is divided into the following 4 areas:

- Patient area
- Quick menu
- Function space and
- Setup

In addition, you will find an information area (information bar or information panel) in the top right-hand corner.

7.1.1 Touchscreen

Operation is designed for use of the touch display, an on-screen keyboard is available to input data.

7.1.2 Camera head buttons

All compatible camera heads and all X-Line videoscopes can be used to call up the programmed functions and for menu control.

- ② Scroll up: Scrolls up the menus, options.
- ③ Menu/selection: Access to the quick menu, with this button (as when pressing the ENTER button) you also activate a marked option on the computer keyboard.
- ④ Scroll down: Scrolls down the menus, options.

7.1.3 Keyboard and mouse

The keyboard and mouse can be used to navigate the menus. A keyboard can be used to enter data.

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
7.2.1 Patient area

Depending on the configuration of the system, the patient area allows one or several patient data sets to be entered manually or for patient data sets to be called up and selected from a patient worklist.

i NOTE: In order to successfully create a patient, the field 'Last name' and (in the case of network storage) the 'Date of birth' must be filled in.

After treatment, images and videos of the current treatment can be viewed in the patient area. More information can be found in chapter 7.7 – Viewing the image material.

Image viewing

Displays a captured still image.
The image does not correspond to the live image.

Video viewing

Shows a recorded video.
The video does not correspond to the live image.

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7.2.2 Quick menu

The quick menu allows you to choose between different camera functions and change the settings of the device. Functions which are not available are shown in gray.

If you operate the device using the camera head buttons or the keyboard, you access the quick menu by pressing the enter button on the keyboard or by pressing the menu button on the camera head/videoscope.

i NOTE: You can change the configuration of the icons in the quick menu via the setup menu (see chapter 7.3).

i NOTE: Some icons in the quick menu are shown in gray until a camera head is connected.

Description of functions of the symbols used

	Full screen The endoscopy image can be shown in full screen mode here (the user interface is no longer shown).
	White balance A white balance is performed here.
i	NOTE: Perform a white balance before every use.
	Freeze You can freeze the image here; during this time the live image is shown in the top right-hand corner of the monitor.
	Zoom This enables digital magnification of the display.
	Alignment The displayed image can be flipped vertically or horizontally and rotated 180° here.

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Enhance

The digital fiberscope filter (anti-moiré/anti-grid filter) can be set to two levels (A, B).

With these filters, interfering artifacts, such as the moiré effect and display of the fiberscope grid, are reduced.

If the honeycombed grid structures are relatively small in the image, Filter A is recommended.

Filter B is particularly recommended in the case of larger honeycomb structures.

Light

The light source can be switched on and off here.

NOTE: If a videoscope with an integrated light source is used, the integrated LED of the videoscope can be switched on and off via this function. In this case, the integrated light source in the TELE PACK + is switched off.

Light source

The light source settings can be displayed and changed here. The light source can be switched on and off, and the mode and light intensity can be set.

Exit

You can exit the quick menu here. Alternatively, you can also exit the menu by clicking the 'Esc' button on the keyboard.

NOTE: The 'Exit' button only appears when using a keyboard or with operation via camera head buttons.

NOTE: Upon pressing the 'Exit' button, the system returns to its previous state.

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Image capture

Still images can be captured here.

NOTE: In order to save the images, a USB storage medium or a configured network storage location are recommended. At the same time, the images can be printed using a connected printer.

NOTE: In addition to or in place of a USB storage medium or network storage, the AIDA® documentation unit can be used to capture still images. These are saved on the AIDA® following the image capture.

Video recording

Video recordings can be started and stopped here.

Training mode

A circle can be displayed in the center of the image here. The training mode can be used specifically for endoscopy training with the TELE PACK +.

Camera brightness

The brightness of the camera can be set here.

Print

Print now: Immediately prints all images in the print queue.

NOTE: This button only appears if a printer is connected.

Swap cameras

Here you can switch between two connected cameras without having to disconnect one camera.

NOTE: This button only appears if two cameras are connected at the same time.

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7.2.3 Function space

Just like the quick menu, the function space (see chapter 7.2) allows you to choose between different camera functions and change the settings of the device (via the setup menu). It automatically contains all functions which were not configured in the quick menu.

i NOTE: Some icons in the function space are shown in gray until a camera head is connected.

The function space is called up using the button shown.

In the setup menu (see 7.2.4) you can access the function space by pressing the setup button (below).

7.2.4 Setup menu

You can change the general TELE PACK+ settings in the setup menu.

Display language
The language of the user interface can be selected here.

Keyboard language
The keyboard language of the on-screen keyboard can be selected here.

i NOTE: The selected keyboard language is also used for a connected, external keyboard.

Stroboscopy settings
The speed of the visible vocal cord movement can be set between 0.2 and 2.5 Hz here. The setting is displayed in Δ Hz when the slow motion mode is activated.

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 **DVI signal output**
You can select the image refresh rate between 50 Hz and 60 Hz for the DVI output here.
 **NOTE:** After changing the frequency, the device must be restarted.

 **Date & time**
The date, time and display format can be set here.

 **Patient management**
The settings for handling patient data can be managed here. Please see chapter 7.5 for details.

 **Module information**
All the necessary information on the system and the connected camera heads and videoscopes can be called up here.

 **Quick menu settings**
The quick menu can be configured here (see chapter 7.3).

 **Network settings**
The settings for the TELE PACK + can be changed here and the user-defined storage location can be specified.
The network function can be used in two different ways. On the one hand in the 'FileShare' mode. In this mode, the generated data is automatically moved to an FTP server (FTP = File Transfer Protocol). To this end, you will need to configure the 'Local' and 'FileShare' tabs. The 'Worklist' mode is also possible, here you call up and select the patient data from the worklist server and then store it on the worklist server again. Please also note chapter 7.8 for details on network storage functionality.

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 **SW update**
The device software and the connected camera heads and videoscopes can be updated here.

 **Software licenses**
The device contains open source software. A list of the open source software included in the product, and the open source software license terms can be viewed here.

 **Logged events**
System logs
Here you can display the system log which the device records on network transmission and export it either via USB or the selected FTP server.
Audit logs
You can view the audit log here or export it via either USB or the selected FTP server.

 **Factory settings**
The device can be reset to these factory settings here, note the information in chapter 7.4.

 **Access & protection**
You can manage access to sensitive areas of the system here.
You will find more details in chapter 7.4 Access and security concept

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7.2.5 Information area

The information area in the top right-hand corner of the screen comprises the information bar and the information panel (see chapter 7.2).

The information bar offers a continuous overview of the system status. You will find more information on the system status and the connected accessories in the information panel. The information panel is folded down by pressing the 'i' symbol in the system bar (see figure).

Symbols in the information area

	Network connection active
	Network connection disconnected
	Camera
	Footswitch
	Printer
	USB storage
	Internal storage

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Overview of connected camera head/videoscope (incl. button assignment)

If a compatible camera head is connected to the TELE PACK +, you can call up the article number and name as well as the camera head button assignment here.

The camera head button assignment is shown true to original; the assignment can be altered using the 'Settings' button.

i NOTE: The camera head button assignment cannot be changed on all videoscopes, note whether a fixed assignment is already printed on the videoscope.

Overview of a connected footswitch

If a compatible footswitch is connected to the TELE PACK +, you can call up the article number and name as well as the button assignment here.

Overview of storage medium used

If a storage medium is selected and connected to the TELE PACK +, you have the option here of calling up the currently used and free storage space.

Should you use the internal storage, you can transfer the recorded data to a USB storage medium here.

i NOTE: Internally saved data is deleted after successful transfer to the USB storage medium in order to create space for new recordings.

You will find more information on the device storage options in chapter 7.4.3 Storage functionality

Overview of a connected printer

If a compatible USB printer is connected to the TELE PACK +, you have the option here of calling up the printer configuration and changing this if necessary. The images per page and the number of copies can be changed and it can also be set whether images are to be printed out automatically during an ongoing treatment.

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The configuration menu is divided into 3 columns. Here functions can be moved, replaced, added or removed from the quick menu.

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- i NOTE:** You can add up to 8 functions to the quick menu.
- i NOTE:** The configuration is saved by restarting the device.

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- ① **Quick menu:**
The current functions are listed here (see 7.2.2).
-  Press this button to add a new function to the quick menu.
- ② **Command:**
The functions in the quick menu can be configured using the 'Command' buttons. To this end, select a function in the (left) column ①.
-  Press this button to move the selected function upwards in the quick menu.
-  Press the 'Delete' button to remove the selected function from the quick menu.
-  Press the 'Replace' button to select another function (see ③) in place of the current one.
-  Press this button to move the selected function downwards in the quick menu.
- ③ **Select Function:** The functions which are not currently configured in the quick menu are listed here.

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7.4 Access and security concept

The TELE PACK + offers two options for protecting the device and sensitive menu areas from access from unauthorized parties. One of the two options must be selected when configuring the system in order to work with sensitive information such as patient data.

i NOTE: Please note that for security reasons, the device must be reset to the factory settings if the password is lost. As such, in addition to all settings, the patient data is also deleted from the device and the device is reset to the default patient privacy settings.

7.4.1 Encryption of sensitive menus via password

By selecting this option, certain menus are protected by a password. The standard password is [90290245] and must be changed.

7.4.2 Role-based access system via individual user accounts

Upon selecting this option, sensitive areas of the menu as well as device functions can only be selected by certain users. Authorization to call up these menus and functions is linked to the specific roles of the user accounts.

Configuration of the role-based access system

When you configure the role-based access system, you must firstly create an administrator account. Moreover, in the first step you must specify password rules which apply to all further passwords. Then, with the help of the standard password, you can create a new password and subsequently make the further settings.

i NOTE: The 'Administrator' role has full access to all device functions and settings.

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- i NOTE:** In addition to the Administrator, the roles of 'OR Personal' and 'Technician' can also be selected. The 'OR Personal' role enables access to patient data, but not to the patient data security and memory areas. The role of 'Technician' enables access to device settings in order, for example, to perform updates without providing access to patient data or other sensitive areas.
- i NOTE:** All safety-relevant actions are logged and recorded by the system. KARL STORZ assumes no legal liability or responsibility for the loss, accuracy, damage or disclosure of patient data which was managed using TELE PACK +.
- Configuration of further user accounts**
You can create further user accounts, edit existing accounts and also delete accounts under this menu item.
A name, corresponding role and a password must be specified for each user account.
 - i NOTE:** Limitation of the number of user accounts.
 - Login message**
You can enter a customized message here which is shown in the login area before granting access to sensitive areas.
 - Password rules**
Here you can set whether and which rules are to be applied when assigning passwords.
 - Automatic logout**
You can set here whether the system should log out a user after expiry of an individually set period of inactivity.
 - Deactivation**
Upon deactivating the role-based access system, the encryption of sensitive menus via password is automatically activated, i.e., an individual password must be created.

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Use of the role-based access system

After you have activated the role-based access system, a further menu item with the name 'Login' appears in the main screen next to the information area.

Here all system users can log in using their user account and password. If no users are logged in, certain areas such as the patient area and the settings for patient data security and memory cannot be called up.

7.4.3 Memory functionality

The TELE PACK + comprises 3 storage options:

- Internal storage
- USB storage and
- Network storage

Depending on the configuration of the recording and storage options in the setup menu under the item Access and security, images and videos with patient data can be captured in the folder and file name. The type of storage can also be selected here.

For storage on USB and the network, the internal storage is used for buffering. This means data is firstly buffered in the internal storage and then transferred to the USB or network directory.

i NOTE: The transfer status of the data can be checked in the information area, see chapter 7.2.5 Symbols in the information area.

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In order to store recordings without any loss of data, the 'Complete treatment' button in the patient area must be pressed. The pending data is then transferred and saved properly. If no network storage location has been configured or no USB storage medium is available, the internal storage of the TELE PACK + can also be used. Approx. 50 GB of storage space are available for internal storage.

i NOTE: Internally stored recordings can only be removed from the internal storage via USB export. Please note that if large volumes of data are saved in the internal storage, a correspondingly large USB storage medium must be available.

i NOTE: The transfer of large data volumes to a USB storage medium can take some time. The current status of data transfer can be checked in the information area.

7.5 Patient data security

i NOTE: The TELE PACK + system does not contain patient data as standard. The patient data administration functions (e.g., entry, access, display, storage and transfer) must be activated by the user.

i NOTE: Upon activating the patient data administration functions, the user assumes responsibility for all patient-related data which are managed using TELE PACK + as per the applicable national data protection laws.

i NOTE: KARL STORZ assumes no legal liability or responsibility for the loss, accuracy, damage or disclosure of patient data which is managed using TELE PACK +.

In the Patient Data Privacy settings you can configure all settings related to the patient information.

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Patient information of any kind is only available if the corresponding settings have been made.

The configuration of the 'Patient Data Privacy' menu is password-protected as standard. You will find further information on password protection in chapter 7.4. Access and security concept.

The 'Patient Data Privacy' menu is divided into 3 areas:

- Patient handling
- Labeling and
- Export mode

7.5.1 Patient handling

1) Worklist setup

- No Worklist / Patient Entry (selected as standard)
- Local Worklist
- Worklist Provider

Here you can establish whether patient data may be used (this is not permitted as standard) and how patient data can be processed with the TELE PACK +. Either via manual entry or via a patient worklist using a network worklist provider (see chapter 7.8 Network communication).

i NOTE: When you are working without entering patient data, the data in the info panel must be exported via 'Transfer to USB'.

2) Capture tools

- No Capture
- Capture Allowed (selected as standard)

Here you can specify whether and how images and videos are to be captured. The following options are available: **no capture, capture of images and videos labeled with patient information and the anonymous capture of media files.**

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3) Storage limit

- Infinite
- Power cycle
- 1 day
- 7 days (selected as standard)
- 30 days
- 90 days

Here you can specify how long data sets are to be stored on the device.

i NOTE: Please note that both patient data sets in the worklist and media files (recorded images and videos) are deleted if this is set here. The data can no longer be restored.

7.5.2 Labeling

The labeling specifies whether the digital recordings on the USB stick or FTP server or the print-outs are to be labeled with patient information or not.

Labeling of the digital recordings:

- No patient data in the patient folder (selected by default)
- Patient data in the patient folder

Labeling of the print-outs:

- No patient data on print-outs (selected by default)
- Patient data on print-outs

7.5.3 Export mode
Selecting the Export mode

- No export
- USB (or USB anonymous) export (selected as standard)
- Network export
- Network and USB export

Here you can specify what type of data export is allowed. In addition to the options of saving data on a USB storage medium or transferring it to a network storage path, it is also possible to enable both storage options.

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i NOTE: The stroboscopy information is displayed once the stroboscopy footswitch is actuated.

i NOTE: The microphone set (40160031) displays the volume in dB, in accordance with airborne sound measurements. This value may deviate from the structure-borne sound measurements.

7.6.1 Connecting a stroboscopy microphone

Insert the stroboscopy microphone 20 1400 30 in the side stroboscopy microphone input socket.

7.6.2 Activating the stroboscopy functionality

The stroboscopy function is activated by connecting the stroboscopy footswitch 40160033 to the USB port.

7.6.3 Starting/stopping the stroboscopy

When the footswitch is operated in the presence of a microphone signal, stroboscopy mode is automatically started and the stroboscopy specific parameters are displayed. Releasing the stroboscopy footswitch ends the stroboscopy mode and reactivates the continuous light mode.

7.6.4 Still image function

The still image function is active in the upper area of the pedal movement range, meaning that the vibration of the vocal cords appears as a static image.

The respective phase or the position of the vocal cords can be adjusted within a range using the footswitch from 0 to 360 degrees.

7.6.5 Slow motion function

If the footswitch is pressed further down beyond the pressure point, slow motion function is activated.

The vocal cords appear in the preset slow motion frequency. The frequency can be set from 0.2 to 2.5 Hz under Settings – Stroboscopy – Slow motion.

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7.6.6 Starting/stopping a recording

The video recording can be started or stopped by pressing the stroboscopy footswitch twice in quick succession.

7.7 Viewing the image material (as of SW version 2.0)

After the treatment, recorded images and videos can be viewed before the current treatment is completed. An overview of the recorded images and videos can be found in the patient area. Images and videos can be viewed on the full screen, selected for printing or deleted.

When viewing the videos, in addition to playing the videos with sound, it is also possible to take snapshots from the video or to view the frames of the video individually. Each scene of the video can be selected using the slider to ensure simple operation.

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7.8 Network communication

The TELE PACK + can be operated in network mode. This operating mode enables the recorded data to be stored at a user-defined storage location when a network connection is active. Two storage locations are available. You can store the data on a standard FTP server or on a worklist server.

Data (images and videos) generated during the examination is temporarily stored in the internal memory with encrypted assignment, and is promptly transferred to the defined server when the network is available.

If the network is not available, all data will be stored internally for the time being. The TELE PACK + sends the data automatically to the selected server as soon as the network becomes available again.

7.8.1 File share mode

In File Share mode a new patient data set is created in the patient area or selected from an existing list. All video data and images generated during the treatment are automatically transferred to the server directory you defined (see chapter 7.9.3). The TELE PACK + will create a subdirectory for the patient.

As soon as a file has been successfully transferred from the internal memory of TELE PACK + to the server directory and treatment is complete, the local copy is deleted.

If the network connection is interrupted during treatment, the TELE PACK + will automatically transfer all data as soon as the network connection becomes available again.

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
7.9 Setting up the network connections

7.9.1 The network dialog

You can call up the network dialog of the TELE PACK + via the setup menu --> Network settings. You can change the settings for the TELE PACK + and specify the user-defined storage location in the network dialog.

The network function can be used in two different ways. One of these is the 'File Share' mode. In this mode, the generated data is automatically moved to an FTP server. To this end, you need to configure the 'LAN' and 'FTP/FileShare' tabs. The other way of using the network function is the Worklist mode, which allows you to access and select the patient data from the worklist server, and to store it on a worklist server again. As two different worklist servers can be used here for access and storage, you will need to have configured all three tabs. Details on this can be found in chapter 7.8 'Network communication'.

i NOTE: In order to use the Worklist functionality, a connection to the KARL STORZ server software SCENARA®.connect is necessary. For further information, please contact your responsible KARL STORZ employee on site. *StreamConnect for US customers

7.9.2 Local settings

The network configuration of the TELE PACK + is defined in the 'LAN' tab. This is performed either automatically via DHCP or manually. If you would like to use the automatic address assignment, you will need to activate DHCP.

If you do not want to automatically configure the TELE PACK +, you will need to enter the necessary data manually.

The configuration set will then need to be tested by pressing the 'Connect' button. The displayed test report confirms successful configuration. The network function is now available.

i NOTE: The TELE PACK + only supports IP version 4.

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
7.9.3 File Share settings

In order to define the storage location for your data in the network, you will need to save the configuration for the server both in the File Share mode and in the Worklist mode in the 'FTP/ FileShare' tab (see figure on the left).

After entering the data, the set configuration must be tested by pressing the 'Connect' button. The displayed test report confirms successful configuration. The network function is now available.

File Share	The name of the directory on the FTP server.
Host name/ IP address	IP address or host name of the FTP server.
Port	The port number used to reach the FTP server. Usually 21.
User	If your FTP server requires access, you will need to store the access data.
Password	
Timeout	Maximum period of time that the TELE PACK + tries to establish a connection to the FTP server.

i NOTE: It is recommended that a user be created for the TELE PACK + on the FTP server.

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
7.9.4 Worklist settings

If you are using a worklist server as a user-defined storage location, you will need to save the configuration for the server in the 'Worklist' tab. The tab 'FTP/FileShare' must also be configured in this case.

After entering the data, the set configuration must be tested by pressing the 'Connect' button. The displayed test report confirms successful configuration. The network function is now activated in the 'Worklist' mode.

Max worklist items count:	Maximum number of entries received by the TELE PACK +.
----------------------------------	--

i NOTE: If this value is set too high, this can have a noticeable impact on the performance of the TELE PACK +. A maximum value of 1000 is recommended here.

Device identity	A unique name for the TELE PACK +
Host name/ IP address	IP address or host name of the worklist server.
Port	The port number used to reach the worklist server is usually 443.
User	If your worklist server requires access, you will need to store the access data.
Password	
Timeout	Maximum period of time that the TELE PACK + tries to establish a connection to the worklist server.

i NOTE: It is recommended that a user be created for the TELE PACK + on the worklist server.

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8.1 Reprocessing

8.1.1 General safety instructions

WARNING: Always unplug the device and disconnect the endoscope before carrying out any maintenance or cleaning.

WARNING: No service and maintenance work may be performed whilst the device is connected to the patient.

CAUTION: Make absolutely sure that no liquid can get into the device housing.

WARNING: Risk of infection: Incorrectly reprocessed medical devices expose patients, users and third parties to a risk of infection as well as the risk that the medical device may malfunction. Observe the manual 'Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments' and the accompanying documentation.

WARNING: When preparing and using the solutions, follow the chemical manufacturer's specifications, paying close attention to proper concentration, exposure time and service life. Prolonged immersion and incorrect concentration may result in damage. Bear in mind the microbiological range of action of the chemicals used.

WARNING: When carrying out any work on contaminated medical devices, the guidelines of the Employers' Liability Insurance Association and equivalent organizations striving to ensure personal safety must be observed.

WARNING: National laws and regulations must be observed.

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
**8.1.2 Wipe-down disinfection of
devices**

Wipe-clean the exterior surfaces of the medical device with a disposable cloth moistened with disinfectant or a ready-to-use soaked disinfectant cloth. Due to their protein-fixating effect and possible material incompatibility, alcohol-based agents must not be used. The chemical manufacturer's specifications regarding material compatibility must be observed. At the end of the necessary exposure time for the disinfectant as specified by the manufacturer, wipe the surface with a dry low lint cloth.

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8.2 Maintenance and safety check

8.2.1 Maintenance

Preventive maintenance is not essential. Regular maintenance can, however, contribute to identifying potential problems before they become serious, thus enhancing the device's reliability and extending its useful service life. Maintenance services can be obtained from your local representative or from the manufacturer.

8.2.2 Safety check

WARNING: *Regardless of the accident prevention regulations or testing intervals for medical devices prescribed in different countries, with this device safety checks/repeat inspections as defined by IEC 62353 must be performed and recorded by a qualified electrician at least once a year.*

Visual inspection

- Check the device and accessories for any mechanical damage which may impair functionality.
- Check that the inscriptions relevant to safety are legible.

Electric measurements

- Protective grounding measured according to IEC 62353 (current edition), and checked for compliance with the limits specified in the standard.
- Touch Current, Earth Leakage Current, and Patient Leakage currents measured according to IEC 62353 (current edition), and checked for compliance with the limits specified in the standard.

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Functional test

- A functional check must be carried out in accordance with the instruction manual (see chapter 6.3.11).
- The safety check must be entered in the unit log-book and the results of the check must be documented.
- If the device does not function reliably and/or safely, it must be repaired.

WARNING: If any defects or shortcomings which could endanger patients, personnel or others are ascertained during these safety checks/repeat inspections, the device must not be operated until such time as these defects or shortcomings have been eliminated by qualified technical servicing.

8.3 Servicing and repair

Defective items of equipment must be serviced and repaired exclusively by persons authorized by KARL STORZ; all repair work must employ original KARL STORZ parts only.

8.4 Disposal

This device has been marked in accordance with the European Directive on waste electrical and electronic equipment (WEEE).

At the end of its useful operating life, dispose of the device as electronic scrap. Please ask KARL STORZ SE & Co. KG, a KARL STORZ subsidiary, or your authorized dealer about the collection point in your area.

Within the scope of application of this Directive, KARL STORZ SE & Co. KG is responsible for the proper disposal of this device.

Control of pollution caused by electronic information products (China RoHS).

NOTE: Please observe the warnings and notes on data security in chapter 7.4 and 7.5 before disposing of the device!

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78579 Neuhausen

Servicehotline: +49 7461 708-980

E-Mail: technicalsupport@karlstorz.com

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8.5 Repairing devices

Only personnel authorized by the manufacturer may repair devices. Interventions described in these instructions for use are exempt from this rule.

In Germany, you can contact the following address for repairs:

KARL STORZ SE & Co. KG

Abt. Reparaturservice

Take-Off Gewerbepark 83

78579 Neuhausen, Germany

Service hotline: +49 7461 708-980

E-mail: technicalsupport@karlstorz.com

In other countries, the respective KARL STORZ branches or specialist dealers are responsible.

8.6 Infection prevention

Contaminated devices may not be shipped. To prevent contact infections and airborne infections, devices must first be decontaminated. KARL STORZ reserves the right to send back contaminated devices.

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8.7 Limitation of liability

KARL STORZ shall be liable for the safe operation, operational reliability and performance of this equipment only subject to the following conditions:

- all assembly operations, system expansions, readjustments, modifications or repairs have been performed by a person or persons duly authorized by KARL STORZ,
- all electrical installations at the location in which the device is connected and operated meet the applicable laws and standards, and
- the device has been used in accordance with its operating instructions at all times.

8.8 Warranty

The warranties provided can be found in the Standard Conditions of Business of KARL STORZ. The medical device must always be sent to your local subsidiary (see chapter 'Subsidiaries'), even during the warranty period.

Opening of the equipment or performance of any repairs or modifications to the equipment by unauthorized persons shall relieve us of any liability for its performance. Any such opening, repair or modification performed during the guarantee period shall void all warranty.

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9 Technical description

9.1 Error messages

Malfunction	Touchscreen warning
LED temperature above 80°C (Message disappears when temperature <75°C again)	Overheating of lamp. Increased lamp temperature.
LED temperature above 90°C (Message disappears when temperature <80°C again)	Overheating of lamp. The lamp will be switched off in a few minutes.
LED temperature above 100°C and lamp switched off by the hardware (Message disappears when temperature <75°C again)	Lamp overheated. Wait a few minutes to allow the device to cool down.
The device does not start up correctly as the fan doesn't work.	The system is switched off due to a fan malfunction – please contact the KARL STORZ Service for further assistance.
LED is switched off due to excess current. The user can attempt to switch the LED on again manually.	The light source was deactivated for safety reasons.
If the temperature in the device increases above 80°C, the message appears, if the temperature goes above 90°C, the device is switched off. If the temperature drops below 75°C, the message is deleted.	Overheating of device. The device will be switched off in a few minutes. Please consult KARL STORZ Service.
The software update could not be carried out for certain reasons.	Upgrade failed. Repeat the procedure after restarting. If the problem persists, please contact a representative or the Service Center.
The camera head is not compatible, the message is shown when changing the camera.	Incompatible head

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9.2 Technical data

TELE PACK +	TP101
Line voltage	100...240 VAC
Power frequency	50/60 Hz
Power consumption	140 VA
Dimensions	470 mm x 387 mm x 222 mm (W x H x D)
Weight	9 kg
Storage/transport temperature	-10°C to +60°C
Humidity (non-condensing)	Operation: 15% ... 80% Storage/transport: 5% ... 85%
Air pressure	Operation: 700 hPa – 1080 hPa Storage/transport: 500 hPa – 1080 hPa
Operating temperature	0°C to 40°C
Operating rated altitude NN	3000 m
Light source	
Lamp	LED
Color temperature	5700 K
Display	
Screen size	18.5"
Resolution	1920 x 1080p
Contrast	1200 : 1
Speaker	
Power	3.2 W

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- Type of protection against electric shock: Protection Class I
- Protection of camera applied part against electric shock:
Depending on the camera connection used, type BF or CF, defibrillation-proof
- Protection of light applied part against electric shock: Type CF, defibrillation-proof
- Protection of microphone for stroboscopy against electric shock: Type BF
- Type of protection against moisture: IPX0

i NOTE: Please observe the electromagnetic compatibility information in chapter 11.

**9.4 Directive compliance
(for TP101)**

According to Medical Device Directive (MDD):

This medical device belongs to Class I
This medical device bears a CE mark in accordance with MDD 93/42/EEC.

**9.5 Software ownership and
licensing**

With respect to KARL STORZ products containing software components, the customer has a nonexclusive, limited, non-transferable license to use the programmed logic, computer programs and/or software supplied by KARL STORZ, in connection with, and incorporated into, KARL STORZ products (Software). The user right is only valid for the delivered form and only for use in accordance with the printed instruction manual from KARL STORZ for the products sold to the customer (and for no other product). The software, and all modifications, enhancements and upgrades thereto, will, at all times, remain the

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property of KARL STORZ. The customer may not duplicate, copy, reverse-engineer, de-compile, or disassemble the software or in any way modify the software. The customer has no right to create derivatives of the software, and the customer may not copy, decipher or re-create the source code of the software. Any and all such modifications or enhancements to the software by the customer, in contravention of this license, will immediately become the sole property of KARL STORZ

The customer hereby acknowledges and agrees (I) that the purchase of KARL STORZ products does not constitute a sale of the software, (II) that the software is the property of KARL STORZ, (III) that the customer neither owns nor acquires any claim or intellectual property right in or to the Software, and (IV) KARL STORZ retains all right, title, and interest in and to the software, at all times, regardless of the form or media in or on which the original or other copies of the software may exist. In the event of a failure of the customer, or its agents, employees or representatives, to comply with any terms and conditions of the License herein granted, the license will, without any further action by KARL STORZ or any other party, immediately terminate.

Attribution notice

This device contains open source software. You can view a list of the open source software included in the product, and the open source software license terms, on the TELE PACK + device itself. Select the tab Software licenses in the setup menu.

To obtain a copy of the source code being made publicly available by KARL STORZ Imaging, Inc. ("KARL STORZ") related to open source software used in this KARL STORZ product, you should send your request, along with the device and serial number, in writing by email to opensource@karlstorz.com or by postal mail to:

KARL STORZ Imaging, Inc.
Attention: Open Source Software Inquiries
1 S. Los Carneros
Goleta, CA 93117
USA

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**10 Spare parts,
recommended accessories**

i NOTE: For use with the product, KARL STORZ recommends the accessories listed in this chapter. If third-party accessories are used, the user must check the correct operation himself.

Not all accessories are offered in all regions.

10.1 Compatible camera heads for X-Line connection

Article/model	Description
TH110	IMAGE1 S™ HX
TH111	IMAGE1 S™ HX-P
TH115**	IMAGE1 S™ D1
TH116	IMAGE1 S™ D1

10.2 Compatible camera heads for C-Line connection

Article/model	Description
TH130**	H1

10.3 Compatible videoscopes for X-Line connection

Article/model	Description
	ENT
11101VP/VN	CCD Video Rhino-Laryngoscope
11101VPS/VNS	Strobo Video Rhino-Laryngoscope
11101HD	HD Video Rhino-Laryngoscope

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Article/model	Description
UROLOGY	
11278VSU/VSUE	Video Uretero-Renoscope FLEX-XC
11278VS/VSE	Video Uretero-Renoscope FLEX-XC
11272VP/VN	Flexible Video Urethro-Cystoscope
11272VPU/VNU	Flexible Video Urethro-Cystoscope
11272VH/VHU	Flexible Video Cystoscope HD-VIEW
BRONCHO	
11900BP/BN*	Video Bronchoscope
11910T/D**	Video Bronchoscope HD
11910P/S**	Video Bronchoscope
10973HD	Video Mediastinoscope
SURGERY	
11292VP/VPU**	Flexible Video Choledochoscope
11292VS/VSU	Flexible Video Choledochoscope (new sensor)
VETERINARY	
60278VS/VSU**	FLEX-XC Veterinary

* needs light adaptor 487VE, adaptor for connecting KARL STORZ videoscopes

10.4 Compatible videoscopes for C-Line connection

Article/model	Description
ANESTHESIA	
8403XS	C-MAC® S Imager
8403XSI	C-MAC® S Imager
8403XSP	C-MAC® S Pediatric Imager

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10.6 Recommended accessories

Article/model	Description
TC001	Videooscope adaptor 0°
TC013**	Videooscope adaptor 90°
TC009	USB to ACC adaptor
TP002	Camera head holder, irrigation bottle holder
20221070	ACC connecting cable, length 180 cm
400F**	Power cord, length 300 cm
400B	'Hospital Grade' power cable (USA)
TC010	Dual-pedal USB footswitch
20040240XX**	USB silicone keyboard with touchpad
20040282	USB stick, 32 GB
20040086	DVI connecting cable, length 200 cm
20040089	DVI connecting cable, length 300 cm
549M**	USB color printer
W21341**	Headset

**Device not currently available in the U.S.

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**10.7 Recommended accessories
light source**

Article/model	Description
487VE	Adaptor for connecting KARL STORZ videoscopes
TL005	Three-way adaptor for light source

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**10.8 Recommended accessories
stroboscopy**

Article/model	Description
40160033	Stroboscopy footswitch
40160031	Microphone set, consisting of:
20140030	Microphone
20140030MH	Microphone holder
20140030T	Adaptor for structure-borne sound
20140030MB**	Membrane

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11 Electromagnetic Compatibility (EMC) information

The TELE PACK+ TP101 is suitable for use in professional healthcare facilities. Professional healthcare facilities include physician offices, dental offices, limited care facilities, freestanding surgical centers, freestanding birth centers, multiple treatment facilities, hospitals (emergency rooms, patient rooms, intensive care, surgical rooms, outside the high-energy shielded room of an ME system for magnetic resonance imaging).

WARNING: Use of this equipment adjacent to or stacked with other equipment should be avoided as it could result in improper operation. If such use is necessary, this and the other devices should be monitored to ensure that they are functioning correctly.

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Accessories and lines for EMC compliance				
Type	Shielded	Length [m]	Ferrite	Use
Power cord	No	3	No	Power supply
Dual-pedal footswitch	Yes	2	No	Unit control, socket rear
USB keyboard with touchpad	Yes	2	No	Data entry
Microphone	No	2	No	Speech recording
Camera	Yes	3	No	Image transmission
DVI cable	Yes	3	Yes	Signal transmission

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 **WARNING:** The use of an accessory, transducer or cable with the TELE PACK + other than those specified in this manual may result in increased emissions or decreased immunity of the TELE PACK +. When using an accessory, transducer or cable with TELE PACK + other than those specified in this manual, it becomes the responsibility of the user of the TELE PACK + to determine compliance with IEC 60601-1-2 when using this item.

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Electromagnetic Compatibility (EMC) information

WARNING: Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the TELE PACK +, including cables specified by the manufacturer. Otherwise, performance may be impaired.

NOTE: The emission characteristics of this equipment make it suitable for use in industrial areas as well as in hospitals (CISPR 11 Class A) and other professional healthcare environments. If the device is used in a domestic environment (for which CISPR 11 Class B is normally required), the device may not offer sufficient protection for radio transmission operation. The user might need to take mitigation measures, such as relocating or re-orienting equipment.

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Table 1 – Compliance level for immunity tests Guidelines and manufacturer's declaration – electromagnetic immunity			
The TELE PACK + is designed for operation in the electromagnetic environment defined below. The user of the TELE PACK + should make sure that it is used in such an environment.			
Interference immunity tests	EN/IEC 60601 test level	Conformity level	Electromagnetic environment – guidelines
Electrostatic discharge (ESD) acc. to IEC 61000-4-2	± 8 kV contact discharge ± 15 kV air discharge	± 8 kV contact discharge ± 15 kV air discharge	Floors should be made of wood, concrete or covered with ceramic tiles. If floors are covered with synthetic material, the relative humidity must be at least 30%.
Electrical fast transient / burst acc. to IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines 100 kHz repetition	± 2 kV/1 kV ± 1 kV for input and output lines 100 kHz repetition	The power supply quality should be that of a typical commercial or hospital environment.
Surges acc. to IEC 61000-4-5	± 1 kV voltage outer conductor – outer conductor ± 2 kV voltage outer conductor – ground	± 1 kV voltage outer conductor – outer conductor ± 2 kV voltage outer conductor – ground	The power supply quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations of the supply voltage acc. to IEC 61000-4-11	<u>Voltage dip:</u> Dip to 0% for 1 cycle @ 0° phase angle Dip to 70% for 25/30 cycles @ 0° phase angle Dropout to 0% for 0.5 cycle @ 0°, 45°, 90°, 135°, 180°, 225°, 270° & 315° phase angles <u>Voltage interruption:</u> 100% for 250/300 periods	<u>Voltage dip:</u> Dip to 0% for 1 cycle @ 0° phase angle Dip to 70% for 25/30 cycles @ 0° phase angle Dropout to 0% for 0.5 cycle @ 0°, 45°, 90°, 135°, 180°, 225°, 270° & 315° phase angles <u>Voltage interruption:</u> 100% for 250/300 periods	The power supply quality should be that of a typical commercial or hospital environment. If the user of the TELE PACK + requires continued operation during power mains interruptions, it is recommended that the TELE PACK + be powered from an uninterruptible power supply or a battery.
Magnetic field at the supply frequency (50/60 Hz) acc. to IEC 61000-4-8	30 A/m at 50 Hz/60 Hz	30 A/m at 50 Hz/60 Hz	If image distortion occurs, it may be necessary to position the TELE PACK + further from sources of electromagnetic fields or to install magnetic shielding. Before installing the device, the electromagnetic field should be measured to ensure that it is sufficiently low.
Interference to immunity test acc. to IEC 61000-4-3 for high-frequency electromagnetic fields	3 V/m 80 MHz to 2.7 GHz *Refer to table 2 for wireless Proximity RF field test levels	3 V/m 80 MHz to 2.7 GHz	
Immunity to conducted disturbances, induced by high-frequency fields acc. to IEC 61000-4-6	3 Vrms to 150 kHz to 80 MHz 1 kHz 80% AM modulation 6 Vrms in ISM band	3 Vrms to 150 kHz to 80 MHz 1 kHz 80% AM modulation 6 Vrms in ISM band	

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Table 3 – Test levels for Radiated and Conducted Immunity Tests
Guidelines and manufacturer's declaration – electromagnetic immunity

The TELE PACK + is designed for operation in the electromagnetic environment defined below. The user of the TELE PACK + should make sure that it is used in such an environment.

Interference immunity tests	EN/IEC 60601 test level	Conformity level	Electromagnetic environment – guidelines
Conducted RF disturbances acc. to IEC 61000-4-6	3 V _{rms} 150 kHz to 80 MHz	3 V _{rms}	<p>Portable and mobile RF communications equipment should be used no closer to any part of the TELE PACK +, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance:</p> $d = 1.2 \sqrt{P}$ <p>Where P is the rated power of the transmitter in watts [W] according to the information provided by the transmitter manufacturer and d is the recommended separation distance in meters [m].</p> <p>Field strength of fixed RF transmitters, as determined by an electromagnetic site survey ^a, should be less than the compliance level in each frequency range ^b.</p> $d = 1.2 \sqrt{P} \quad 80 \text{ MHz to } 800 \text{ MHz}$ $d = 2.3 \sqrt{P} \quad 800 \text{ MHz to } 2.5 \text{ GHz}$
Radiated RF disturbances acc. to IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	<p>Interferences may occur in the vicinity of equipment marked with the following symbol:</p> 

Note: At 80 MHz and 800 MHz, the higher frequency range applies.

Note: These guidelines may not apply in all situations. The propagation of electromagnetic waves is affected by absorptions and reflections of buildings, objects, and people.

^a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast, and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed transmitters, an electromagnetic site survey should be considered. If the measured field strength at the location where the TELE PACK + is used exceeds the above compliance levels, the TELE PACK + should be monitored to ensure proper function. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the TELE PACK +.

^b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

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Table 4 – Emission class and group		
Guidelines and manufacturer's declaration – electromagnetic emissions		
The TELE PACK + is designed for operation in an environment as defined below. The customer or user of the TELE PACK + should ensure that it is used in such an environment.		
RF emissions acc. to CISPR 11	Group 1	The TELE PACK + uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions acc. to CISPR 11	Class A	The TELE PACK + is suitable for use in all establishments other than domestic and those directly connected to the public low voltage power supply network that supplies buildings used for domestic purposes.
Emission of harmonic oscillations acc. to IEC 61000-3-2	Class A	
Voltage fluctuations/flicker emissions as per IEC 61000-3-3	complies	

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Table 5

Recommended separation distances between portable and mobile RF communications equipment and the TELE PACK +

The TELE PACK + is intended for use in an electromagnetic environment in which the RF disturbances are controlled. The customer or user of the TELE PACK + can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the unit as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter [W]	Separation distance d [m] according to frequency of transmitter		
	150 kHz to 80 MHz $d = 1.2 \sqrt{P}$	80 MHz to 800 MHz $d = 1.2 \sqrt{P}$	800 MHz to 2.5 GHz $d = 2.3 \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters whose maximum rated power is not listed in the table above, the recommended separation distance d in meters (m) can be estimated using the equation from the respective column, whereby P is the maximum rated power of the transmitter in watts (W) according to the transmitter manufacturer.

Note: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

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Compatibility with surgical HF instruments

This device was tested for compatibility with HF surgical devices. It has been tested according to IEC 60601-2-2 Annex BB.

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KARL STORZ Endoscopy Canada Ltd.
7171 Millcreek Drive, Mississauga, Ontario
L5N 3R3 Canada
Phone: +1 905 816-4500, Fax: +1 905 816-4599
Toll free phone: 1-800-268-4880 (Canada only)
Toll free fax: 1-800-482-4198 (Canada only)
E-Mail: info-canada@karlstorz.com

KARL STORZ Endoscopy-America, Inc.
2151 East Grand Avenue
El Segundo, CA 90245-5017, USA
Phone: +1 424 218-8100, Fax: +1 424 218-8525
Toll free phone: 800 421-0837 (USA only)
Toll free fax: 800 321-1304 (USA only)
E-Mail: communications@ksea.com

KARL STORZ Veterinary Endoscopy-America, Inc.
1 South Los Cameros Road
Goleta, CA 93117, USA
Phone: +1 805 968-7776, Fax: +1 805 685-2588
E-Mail: info@karlstorzvet.com

KARL STORZ Endoscopia Latino-America, Inc.
815 N. W. 57th Avenue, Suite 480
Miami, FL 33126-2042, USA
Phone: +1 305 262-8980, Fax: +1 305 262-8986
E-Mail: info@ksela.com

KARL STORZ Endoscopia México S.A. de C.V.
Edificio Atlántic, Oficina 3G
Calle D e/ 1ra y 3ra
10400 Vedado, Havana, Cuba
Phone: +537 836 95 06, Fax: +537 836 97 76
E-Mail: kstorzcuba@gmail.com

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E-Mail: infovietnam@karlstorz.com
KARL STORZ Endoscopy China Ltd.
Room 2503-05, 25F AXA Tower, Landmark East,
No. 100 How Ming Street, Kwun Tong, Kowloon,
Hong Kong, People's Republic of China
Phone: +852 28 65 2411, Fax: +852 28 65 4114
E-Mail: inquiry@karlstorz.com.hk

KARL STORZ Endoscopy (Shanghai) Ltd., Beijing Branch
Room 1805-1807, Building B, 18F Beijing IFC
No. 8, Jianguomenwai Street, Chaoyang District,
100022, Beijing, People's Republic of China
Phone: +86 10 5638188, Fax: +86 10 5638199
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd., Shanghai Branch
Room 701A Building 5 & Room 501 Building 7,
No. 3000 Longdong Avenue, Pilot Free Trade Zone, 201203,
Shanghai, People's Republic of China
Phone: +86 21 60339888, Fax: +86 21 60339808
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd., Chengdu Branch
Room 803-805, 8F Jin Jiang International Building No. 1 West Linjiang Road, Wuhou District,
6100414, Chengdu, People's Republic of China
Phone: +86 28 86587977, Fax: +86 28 86587975
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd., Shenyang Branch
Room 2001-2005, 20F N-MEDIA International Center,
No. 167 Youth Avenue, Shenhe District,
110014, Shenyang, People's Republic of China
Phone: +86 24 23181118, Fax: +86 24 23181119
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd., Guangzhou Branch
Room 02B & 03 & 04A, 35F Teem Tower,
No. 208 Tianhe Road, Tianhe District,
510620, Guangzhou, People's Republic of China
Phone: +86 20 87321281, Fax: +86 20 87321286
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy Asia Marketing Pte Ltd.
No. 8 Commonwealth Lane #03-02
Singapore 149555, Singapore
Phone: +65 69229150, Fax: +65 69229155
E-Mail: infoasia@karlstorz.com

KARL STORZ Endoscopy Singapore Sales Pte Ltd
No. 8 Commonwealth Lane #03-02
Singapore 149555, Singapore
Phone: +65 69229150, Fax: +65 69229155
E-Mail: infoasia@karlstorz.com

KARL STORZ SE & Co. KG
Representative Office Indonesia
Sinarmas MSIG Tower Level 37
Jl. Jend. Sudirman No. Kav. 21
Jakarta Selatan
DKI Jakarta 12920
E-Mail: infoindonesia@karlstorz.com

KARL STORZ Endoscopy Korea Co. Ltd.
9F Hyowon-Building
97, Jungdae-ro, Songpa-gu
05719 Seoul, Korea
Phone: +82-70-4350-7474, Fax: +82-70-8277-3299
E-Mail: infokorea@karlstorz.com

KARL STORZ Endoscopy Taiwan Ltd.
12F, No. 192, Sec. 2, Chung Hsin Rd.,
Sindian District, New Taipei City, Taiwan
Phone: +886 933 014 160, Fax: +886 2 8672 6399
E-Mail: info-tw@karlstorz.com

KARL STORZ SE & Co. KG
Representative Office Philippines
1901 Picadilly Star Bldg., 4th Avenue, BGC
Taguig City 1636, Philippines
Phone: +63 2 317 45 00, Fax: +63 2 317 45 11
E-Mail: philippines@karlstorz.com

KARL STORZ Endoscopy Japan K. K.
Stage Bldg. 8F, 2-7-2 Fujimi
Chiyoda-ku, Tokyo 102-0071, Japan
Phone: +81 3 6380-8622, Fax: +81 3 6380-8633
E-Mail: info-jp@karlstorz.com

KARL STORZ Endoscopy New Zealand Ltd.
31 Morningside Drive Mt Albert
Auckland, 1025, New Zealand
PO Box 56 511, Dominion Rd
Auckland, 1446, New Zealand
Phone: +64 9 846 6044, Fax: +64 9 846 6808
Toll free: +64 508 84 84 84 (New Zealand only)
E-Mail: sales@karlstorz.co.nz

KARL STORZ Endoscopy Australia Pty. Ltd.,
68 Waterloo Road, Macquarie Park NSW 2113, Australia
PO Box 50 Lane Cove NSW 1595, Australia
Phone: +61 (0)2 9490 6700
Toll free: 1800 996 562 (Australia only)
Fax: +61 (0)2 9420 0695
E-Mail: karlstorz@karlstorz.com.au
www.karlstorz.com

STORZ
KARL STORZ—ENDOSKOPE

KARL STORZ SE & Co. KG

Dr.-Karl-Storz-Straße 34
78532 Tuttlingen

Postfach 230
78503 Tuttlingen
Germany

Telefon: +49 7461 708-0
Telefax: +49 7461 708-105
E-Mail: info@karlstorz.com
Web: www.karlstorz.com