# STORZ

## KARL STORZ — ENDOSKOPE

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IMAGE1 S™ 4U Rubina, OPAL1® NIR/ICG Kamerakopf TH121

### INSTRUCTION MANUAL

IMAGE1 S™ 4U Rubina, OPAL1® NIR/ICG camera head TH121

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Thank you for placing your trust in the name KARL STORZ. Like all of our other products, this product is the result of years of experience and great care in manufacture. You and your organization have decided in favor of a modern, high quality piece of equipment from KARL STORZ.

This instruction manual, in conjunction with the instruction manual for the IMAGE1 S™ (96206286D), is intended to serve as an aid in the proper connection and operation of the camera head. All required details and all actions on your part are clearly presented and explained. Please read this instruction manual carefully before proceeding with the equipment. Keep this instruction manual available for ready reference in a convenient location near the equipment.

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**It is recommended to check the suitability of the product for the intended procedure prior to use.**

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## 1 Safety instructions

Safety instructions are measures intended to protect the user and patients from the risks which could arise through the use of the system.

### 1.1 Explanation of warnings and cautions

Please read this manual and follow the instructions carefully. The words **Warning**, **Caution** and **Note** convey special meanings. Wherever they are used in this manual, the accompanying text should be carefully reviewed to ensure the safe and effective operation of the system. To make these words stand out more clearly, they are accompanied by a pictogram.

**WARNING:** A Warning indicates that the personal safety of the patient or physician is at risk.

**CAUTION:** A Caution indicates that particular service procedures or precautions must be followed to avoid possible damage to the device.

**NOTE:** A Note indicates special information about operating the device or clarifies important information.

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**1.2 General safety information**

**WARNING:** Read this instruction manual thoroughly before using the equipment.

**WARNING:** Check that the equipment is functioning properly prior to each surgical intervention. Ensure that the correct video image appears on all video monitors before beginning each intervention. It is recommended that a replacement system is kept available.

**WARNING:** Perceived electric shocks may be experienced when using the camera system in conjunction with electrosurgical units.

**WARNING:** When using high frequency surgical accessories, care must be taken to avert contact between the HF surgical device and the endoscope or endoscopic camera. Inadvertent contact may represent a thermal hazard to the surgeon.

**WARNING:** When energized endoscopes are used with energized endotherapy devices, patient leakage currents may be additive. It should be noted that this is particularly important if a Type CF applied part endoscope is used, in which case a Type CF applied part energized endotherapy device should be used in order to minimize total patient leakage current.

**WARNING:** Always ensure the connector is completely clean and free of debris. If moisture is present, dry thoroughly with sterile towel prior to insertion.

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-  **WARNING:** HF-generating equipment is dangerous and cannot be allowed to make contact with the camera head or system via the accessories or endoscopes connected to the camera head or system.
-  **WARNING:** Risk of blindness. Never look into the open end of a light cable.
-  **WARNING:** Risk of burns. Contact with the open end of a light cable connected to a light source may cause burns.
-  **WARNING:** Test all equipment prior to each surgical use. In the event that the image becomes unusable during surgery, the camera may be disengaged from the endoscope and the procedure continued optically. If this is not possible, it is up to the discretion of the surgeon how best to proceed. Availability of a spare system is recommended.
-  **WARNING:** High energy radiated light through endoscopes may give rise to high temperatures in front of the light outlet and to the tip of the endoscope. To minimize the risk of burns, always adjust the light source to the minimum illumination intensity necessary to achieve optimum illumination of the endoscopic scene when coupled to the video camera.

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## 2 Intended Use Camera Heads

KARL STORZ camera heads, in combination with the corresponding camera control unit, light source, monitor and telescopes or microscopes, are used for real-time display of the endoscopic or open operating field in diagnostic or surgical procedures.

Additionally, they are used for detection during fluorescence imaging in the visible and near-infrared spectral range.

KARL STORZ camera heads have no direct contact with the patient.

### 2.1 Contraindications

No contraindications relating directly to the medical device are currently known. The responsible physician must decide whether the foreseen application is admissible based on the general condition of the patient.

### 2.2 User qualification

The IMAGE1 S™ 4U Rubina, OPAL1® NIR/ICG camera head may only be used by physicians and medical assistants who have a corresponding specialized qualification and who have been instructed in the use of the IMAGE1 S™ 4U Rubina, OPAL1® NIR/ICG camera head.

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


- ① Instrument coupler
- ② Focusing ring/knob
- ③ User control button: for scrolling (downwards) through the menu or activating a camera function
- ④ User control button: for scrolling (upwards) through the menu or activating a camera function
- ⑤ User control button: for accessing the menu/ making a selection

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**4 Symbols employed**

	Fragile, handle with care
	Keep dry
	Temperature limit
	Serial number
	Number of products in the product packaging
	CE marking
	top
	Consult instructions for use
	Manufacturer
	Federal (USA) law restricts this device to sale by or on the order of a physician.

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**5 Operating instructions  
for the IMAGE1 S™ 4U  
Rubina, OPAL1® NIR/ICG  
camera head**

**i** **NOTE:** While operating with surgical lasers intensive laser light may occur and affect image visualization.

**5.1 Frequency of the IMAGE1 S™  
4U Rubina, OPAL1® NIR/ICG  
camera head (50/60 Hz)**

The IMAGE1 S™ 4U Rubina, OPAL1® NIR/ICG camera head is suitable for operation at either 50 Hz or 60 Hz. The camera head must be selected to suit the local conditions.

**i** **NOTE:** The method for switching frequency is described in the instruction manual IMAGE1 S™ – TC200, TC201, TC300, TC301, TC304 (96206286D).

**5.2 Operating the IMAGE1 S™  
4U Rubina, OPAL1® NIR/ICG  
camera head**

**5.2.1 Connecting an endoscope  
and light cable**

The IMAGE1 S™ 4U Rubina, OPAL1® NIR/ICG camera head has an integrated instrument coupler and integrated telescope. To connect an endoscope to the instrument coupler, turn the outer ring of the coupler clockwise and insert the endoscope eyepiece. To prevent rotation of the endoscope, the outer ring of the instrument coupler must be turned counterclockwise until tight. Connect the light cable to the endoscope eyepiece by tightening the knurled screw on the light cable through a quarter turn.

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**CAUTION:** Heat from the light source can cause disinfectant residue on the light cable or on the light inlet of the endoscope to burn. This can have a detrimental effect on the light conduction and thus the image quality.

**NOTE:** For further information regarding light cables, please refer to the manual 'Fiberoptic light cables' (97000213), which can be downloaded or requested at [www.karlstorz.com](http://www.karlstorz.com).

### 5.2.2 Adjusting the focus

The focus on the IMAGE1 S™ 4U Rubina, OPAL1® NIR/ICG camera head can be adjusted by turning the focusing ring/knob ② on the camera lens.

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IMAGE1 S™ 4U Rubina, OPAL1® NIR/ICG camera head	
Image sensor	1/2.5"
Image format	16:9
Image refresh rate	50 Hz/60 Hz
Focal length	19 mm
Dimensions (L x H x W)	150 mm x 55 mm x 41 mm
Weight	260 g

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**6.3 Standard compliance**
**According to IEC 60601-1; EC 60601-2-18.**

Provides Type CF protection when used with KARL STORZ Camera Control Units that have CF symbol on connection port.

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**6.4 Directive compliance**
**According to the Medical Device Directive (MDD):**

This medical device belongs to Class I

This medical device bears a CE mark in accordance with MDD 93/42/EEC.

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**7 Reprocessing  
(not validated for  
US customers)**

In the following chapter, you will find information on care and cleaning of the IMAGE1 S™ 4U Rubina, OPAL1® NIR/ICG camera head.

**7.1 Care and handling**

To maximize the life of your camera, please follow the guidelines given below for the correct handling and care of the camera head:

DO NOT PLUG THE CAMERA HEAD CABLE INTO THE PROCESSOR RECEPTACLE IF IT IS MOIST OR NOT CLEAN. Always ensure that the connector is completely clean. Remove any moisture present with a sterile cloth before connecting.

DO NOT PICK UP THE CAMERA HEAD BY ITS CABLE. Always pick up the camera head first, and then pick up and support the cable and connector. The cable contains many tiny conductors and must therefore be handled with care.

DO NOT STORE THE CAMERA SYSTEM in direct sunlight or at high temperatures.

ENSURE THAT THE CABLE IS NOT COILED TOO TIGHTLY OR FOLDED UP WHILE THE CAMERA IS IN STORAGE. Always arrange the cable in loose coils with a diameter of 6 inches (15 cm) or more.

DO NOT PULL FORCEFULLY ON THE CABLE WHEN UNCOILING. Undesirable 'kinks' may result. Always allow the cable to uncoil gently and loosely.

HANDLE THE SYSTEM WITH CARE, avoiding any strong impacts. When transporting, use the original packaging to prevent damage to the device.

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**Reprocessing  
(not validated for  
US customers)**

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**NOTE:** The manual 'Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments' (art. no. 96216003) can be downloaded or requested by visiting [www.karlstorz.com](http://www.karlstorz.com).

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- Brush: 27652

**7.2.4 Manual wipe-down disinfection**

The exterior surfaces of the medical device must be wiped clean with a disposable cloth moistened with disinfectant or a ready to use soaked disinfectant cloth. Due to their protein-fixating effect and material incompatibility, alcohol-based agents must not be used. The chemical manufacturer's specifications regarding material compatibility must be observed. At the end of the necessary exposure time, wipe the surface with a dry low lint cloth.

**7.2.5 Preparation for cleaning and disinfection**

Heavy soiling, corrosive solutions and pharmaceuticals must be removed from the medical device immediately after use. To this end, preclean the medical device by wiping down and rinsing, for example. As a general rule, KARL STORZ recommends manual precleaning under cold running water.

**7.2.6 Manual cleaning**

The medical device must be completely immersed in a cleaning solution. To ensure bubble-free wetting, the lumina in particular must be filled. After the necessary exposure time, clean the instrument mechanically with the aid of brushes or sponges. Finally, it must be rinsed with cold water to ensure neutralization.

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### 7.2.7 Manual disinfection

The medical device must be completely immersed in a disinfectant solution. To ensure bubble-free wetting, the lumina in particular must be filled. At the end of the necessary exposure time, the medical device must be rinsed several times in order to remove all chemical residues. To this end, water of the best possible quality must be used, taking into account the national regulations. Finally, all of the surfaces, joints, openings, channels and lumina are dried completely with (preferably sterile) compressed air in accordance with the national regulations. The cleaning gun with accessories (Art. no. 27660) is ideal for this purpose.

### 7.2.8 Assembly, inspection and care

The cleaned and disinfected medical device must be visually inspected for cleanliness, completeness, damage and dryness:

- If residues or contamination are still present, the medical device must be manually cleaned and subjected to a full cleaning and disinfection procedure once more.
- Damaged or corroded medical devices must be withdrawn from use.
- Afterwards, a functional check must be carried out.

**i** **NOTE:** During care procedures, use items from the catalog 'Care, Sterilization and Storage Techniques'.

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### 7.2.10 Sterilization

The procedures as well as the process-relevant parameters for the individually validated methods are described in detail in the manual 'Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments' (Art. no. 96216003). The method must be selected taking into account the respective applicable national requirements and in consultation with the device and product manufacturers.

The following sterilization methods have been validated and approved by KARL STORZ for this medical device:

#### Hydrogen peroxide (H<sub>2</sub>O<sub>2</sub>) sterilization – ASP STERRAD®

**WARNING:** Medical devices must be sterilized in dismantled condition.

**WARNING:** Remember that there are restrictions concerning the instruments which can be sterilized in the different STERRAD® sterilization systems in respect of lumen dimensions and material.

**CAUTION:** For detailed information, please consult the user handbook of the respective device.

**NOTE:** Refer to the 'STERRAD® Sterility Guide' to ensure that the medical device concerned can be sterilized using the different STERRAD® devices.

**NOTE:** Sterilization is not possible on surfaces which have been lubricated and oiled.

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- This reprocessing step is part of the validated process.
- This reprocessing step may be carried out additionally in order to aid cleaning. The sterilization method is approved in respect of material compatibility, however must be validated by the operator on site in respect of effectiveness.
- x On the grounds of potential damage, this reprocessing step must not be performed.

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- Bürsten: 27652

### 8.2.2 General warnings

**WARNING:** Risk of infection: Incorrectly reprocessed medical devices expose patients, users and third parties to a risk of infection as well as the risk that the medical device may malfunction. Observe the 'Cleaning, Disinfection, Care, and Sterilization of KARL STORZ Instruments' instructions and the accompanying documentation.

**WARNING:** Risk of infection: These medical devices are not sterile when delivered. The use of non-sterile medical devices poses a risk of infection for patients, users and third parties. Inspect medical devices for visible contamination. Visible contamination is an indication that reprocessing has not been carried out, or has been carried out incorrectly. Reprocess the medical devices before initial use and before and after every subsequent use using validated procedures.

**WARNING:** When carrying out any work on contaminated medical devices, the guidelines of the Employers' Liability Insurance Association and equivalent organizations striving to ensure personal safety must be observed.

**CAUTION:** When preparing and using the solutions, follow the chemical manufacturer's specifications, paying close attention to proper concentration, exposure time and service life. Prolonged immersion and incorrect concentration may result in damage. Bear in mind the microbiological range of action of the chemicals used.

**CAUTION:** National laws and regulations must be observed.

### 8.2.3 Accessories

Accessories required for carrying out reprocessing:

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#### 8.2.4 Bedside Pre-cleaning

1. Immediately after a procedure, at the point of use, wipe the camera head with a soft, lint-free, disposable cloth moistened with water or enzymatic cleaning solution to remove all heavy soiling.
2. Contain and transport the camera head to the decontamination area for cleaning.
3. If immediate cleaning is not feasible, place the camera head in a container and soak in neutral/mild pH enzymatic cleaning solution (e.g. Enzol), diluted with distilled/demineralized water to proper concentration per manufacturer's instructions for no longer than an hour.

**i NOTE:** Complete cleaning of the camera head should be started within two hours of the bedside pre-cleaning.

#### 8.2.5 Manual cleaning

1. Prepare neutral pH enzymatic solution (e.g. Enzol) in a container large enough to fully submerge the camera head. Dilute the solution to proper concentration per manufacturer's instructions.
2. Under running cold tap water:
  - Rinse the camera head to remove all gross debris, paying close attention to the seam beneath the focus ring, for 2 minutes.
  - Brush all camera head surfaces using the soft bristle brush (P/N 27652), paying close attention to the seam beneath the focus ring. Brush for 2 minutes or until there is no visible debris on the brush.
- i NOTE:** If heavy debris collects on the brush while brushing, rinse the brush until there is no visible debris on the brush.
3. Completely immerse the camera head in a neutral pH enzymatic solution. Keep immersed for 5 minutes.

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**Reprocessing  
(validated for  
US customers)**

- !** **CAUTION:** Do not brush the glass window of the camera head. Doing so may cause damage to the glass. Use a soft, lint-free cloth.
4. While immersed, use the soft bristle brush (P/N 27652) to:
    - Brush all camera head surfaces, paying close attention to the seam beneath the focus ring. Brush for 1 minute or until there is no visible debris on the brush.
  - i** **NOTE:** If heavy debris collects on the brush while brushing, rinse the brush until there is no visible debris on the brush.
  5. While immersed, use the soft bristle brush (P/N 27652) to brush the following locations for 30 seconds each:
    - Entire silicone cable
    - Card edge connector
    - Camera head buttons
    - Cam-Lock mechanism, including its crevices within the coupling mechanism
  - i** **NOTE:** If heavy debris collects on the brush while brushing, rinse the brush until there is no visible debris on the brush.
  6. While immersed, use a soft, lint-free cloth to wipe the entire camera head cable.
  7. Remove the camera head from the cleaning solution.
  8. Under running cold tap water, rinse the entire camera head and cable for 3 minutes including:
    - Focus ring
    - Camera head buttons
    - Cam-Lock coupling mechanism
    - Crevices and seams
    - Entire cable
  9. Repeat rinsing (step 8) for a total of 3 times.
  10. Dry the camera head with a soft, lint-free, disposable cloth or compressed air.

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## Reprocessing (validated for US customers)

The following STERRAD® sterilization methods have been validated and approved by KARL STORZ for this medical device:

- STERRAD® 100S
  - STERRAD® NX® standard cycle
  - STERRAD® 100NX® standard cycle
1. Place the camera head in an FDA-cleared STERRAD® compatible instrument tray.
  2. Wrap the sterilization tray with two single layers of FDA-cleared polypropylene sterilization wrap, using a sequential double wrapping technique. If using a sterilization sealed container, refer to the manufacturer's instructions for use on using the correct disposable filters for the indicated sterilization method.
  3. Load the STERRAD® sterilizer, arranging the items such that the hydrogen peroxide plasma can surround them. Do not allow any items to touch the wall of the sterilizer.
  4. Start the STERRAD® sterilizer.

### Hydrogen peroxide (H<sub>2</sub>O<sub>2</sub>) sterilization – STERIS® V-PRO®

The following STERIS® V-PRO® sterilization methods have been validated and approved by KARL STORZ for this medical device:

- V-PRO® 1 'Lumen' Cycle
  - V-PRO® 1 Plus 'Lumen' and 'Non-lumen' Cycle
  - V-PRO® maX 'Lumen' and 'Non-lumen' Cycle
1. Place the camera head in an FDA-cleared V-PRO® compatible instrument tray.
  2. Wrap the sterilization tray with two single layers of FDA-cleared polypropylene sterilization wrap, using a sequential double wrapping technique. If using a sterilization sealed container, refer to the manufacturer's instructions for use on using the correct disposable filters for the indicated sterilization method.
  3. Load the V-PRO® sterilizer, arranging the items such that the vaporized hydrogen peroxide can surround them. Do not allow any items to touch the wall of the sterilizer.
  4. Start the V-PRO® sterilizer.

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**Ethylene oxide (EO) sterilization**

1. Place the camera head in an FDA-cleared sterilization tray.
2. Double wrap the sterilization tray with polypropylene wrap. If using a sterilization sealed container refer to the manufacturer's instructions for use on using the correct disposable filters for the indicated sterilization method.
3. Start sterilization under the following parameters:

**Conditioning parameters**

Temperature: 55 °C (131 °F)  
Humidity: 70 % RH  
Conditioning dwell time: 30 minutes

**Sterilization parameters**

Sterilant: 100 % ethylene oxide  
Temperature: 55 °C (131 °F)  
Gas concentration: 735±30 mg/l  
Exposure time: 180 minutes

**Aeration parameters**

Time: 12 hours, minimum  
Temperature: 55 °C (131 °F)

**Chemical low-temperature sterilization with peracetic acid – STERIS® System 1E®**

1. Place the camera head in one of the following trays:
  - C1200E SS1E General Processing Tray/ Container
  - C1220E SS1E Directed Flow Processing Tray/Container
  - C1140E SS1E Flexible Endoscope Processing Tray/Container
  - C1160E SS1E Universal Flexible Processing (open) Tray
2. Load the tray into the SS1E and start the sterilizer.

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**8.2.9 High-level disinfection**  
The following high-level disinfection methods have been validated and approved by KARL STORZ:**Revital-Ox™ RESERT® (e.g. 2.0% Accelerated Hydrogen Peroxide solution)**

In order to prevent potential damage, do not place camera head in a container with other instruments.

Completely immerse the camera head in the undiluted 2.0% hydrogen peroxide solution for a minimum of 8 minutes at 20 °C (68 °F). It must be ensured that the surfaces are covered and that no air bubbles are present.

To finish, immerse the camera head three times in water. Each immersion rinse should be for a minimum of 1 minute in duration. Discard the water after each rinse, as it will be contaminated with the disinfectant solution. Use fresh water for each rinse. Thorough rinsing of the camera head with water is very important for removing any residual disinfectant.

Dry the camera head and the card-edge connector with a lint-free, sterile cloth. Ensure that the card-edge connector is completely clean and free of moisture before inserting it into the camera processor. Wipe the camera head connector with isopropyl alcohol to remove any residual disinfectant.

**8.2.10 Reprocessing limits**

The end of the product's service life is largely determined by wear, reprocessing methods, the chemicals used and any damage resulting from use.

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## 9 Maintenance

### 9.1 Maintenance

Preventive maintenance is not essential. Regular maintenance can, however, contribute to identifying potential problems before they become serious, thus enhancing the instrument's reliability and extending its useful service life. Maintenance services can be requested from your local representative or from the manufacturer.

### 9.2 Servicing and repair

Defective items of equipment must be serviced and repaired exclusively by KARL STORZ or by persons authorized by KARL STORZ; all repair work must employ original KARL STORZ parts only.

**CAUTION:** Do not attempt to open the camera head or lens. Opening these sealed assemblies affects the camera's soakability and voids all warranties. Do not remove the covers on the camera processor. Warranty is voided unless removed by authorized personnel.

### 9.3 Disposal

This unit has been marked in accordance with the European Directive on waste electrical and electronic equipment (WEEE).

At the end of its useful operating life, dispose of the unit as electronic scrap.

Please ask either KARL STORZ SE & Co. KG, a KARL STORZ subsidiary or your specialist dealer for information on your local collection point.

Within the scope of application of this Directive, KARL STORZ SE & Co. KG is responsible for the proper disposal of this unit.

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**9.4 Repairing devices**

Only personnel authorized by the manufacturer may repair devices. Interventions described in these instructions for use are exempt from this rule.

In Germany, you can contact the following address for repairs:

KARL STORZ SE & Co. KG  
Abt. Reparaturservice  
Take-Off Gewerbepark 83  
78579 Neuhausen

Service hotline: +49 7461 708-980  
E-mail: technicalsupport@karlstorz.com

In other countries, the respective KARL STORZ branches or specialist dealers are responsible.

**9.5 Infection prevention**

Contaminated devices may not be shipped. To prevent contact infections and airborne infections, devices must first be decontaminated. KARL STORZ reserves the right to send back contaminated devices.

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### 9.6 Limitation of liability

KARL STORZ shall be liable for the safe operation, operational reliability, and performance of this equipment only subject to the following conditions:

- all assembly operations, system expansions, readjustments, modifications, or repairs have been performed by a person or persons duly authorized by KARL STORZ;
- all electrical installations at the location of use meet applicable national and local electrical codes, and
- the unit has been used in accordance with its operating instructions at all times.

### 9.7 Warranty

The guarantees provided can be found in the Standard Conditions of Business of KARL STORZ.

The medical device must always be sent to your local subsidiary (see ‘Subsidiaries’ section), even during the warranty period.

Opening the equipment or performance of any repairs or modifications to the equipment by unauthorized persons shall relieve us of any liability for its performance. Any such opening, repair, or modification performed during the warranty period shall void all warranty.

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**10 Compatible  
camera control unit (CCU)**

Item	Cat. No.
IMAGE1 S™ 4U-LINK in combination with IMAGE1 S CONNECT® II	TC304 TC201*

\* IMAGE1 camera control units are available with various configurations and image outputs. These camera control units are compatible with IMAGE1 S™ camera heads. Please check your local region for device configuration availability.

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**11 Spare parts,  
recommended accessories**

Please check your local region for recommended accessory availability.

Article	Order no.
Adaptor, autoclavable, allows telescopes to be changed under sterile conditions	533TVA
Integral rotating beamsplitter, for simultaneous viewing through the telescope and on the monitor	9530BD
Dust cap for camera heads	6349190
Sterile drapes* for KARL STORZ cameras	
Instruction Manual	96286011D
Instruction manual IMAGE1 S™	96206286D

\* For safe and comfortable use of KARL STORZ cameras, we recommend the use of sterile drapes.

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**12 Information on  
electromagnetic  
compatibility (EMC)**

The IMAGE1 S™ 4U Rubina, OPAL1® NIR/ICG camera head TH121 used with the IMAGE1 S™ is suitable for use within close proximity to an active HF electro-surgical device in professional healthcare facility environment. Professional healthcare facility includes physician offices, dental offices, limited care facilities, freestanding surgical centers, freestanding birth centers, multiple treatment facilities, hospitals (emergency rooms, patient rooms, intensive care, surgical rooms, outside the RF shielded room of an ME system for magnetic resonance imaging).

**WARNING:** Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the IMAGE1 S™ 4U Rubina, OPAL1® NIR/ICG camera head, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

**NOTE:** The emission characteristics of this equipment make it suitable for use in industrial areas as well as in hospitals (CISPR 11 Class A) and other professional healthcare environment. If it is used in a residential environment (for which CISPR 11 Class B is normally required) this equipment might not offer adequate protection to radio communication service. The user might need to take mitigation measures, such as relocating or re-orienting equipment.

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<b>Table 1 – Compliance level for Immunity tests</b> <b>Guidance and manufacturer's declaration – electromagnetic immunity</b>			
The IMAGE1 S™ 4U Rubina, OPAL1® NIR/ICG camera head is intended for use in the electromagnetic environment specified below. The user of the IMAGE1 S™ 4U Rubina, OPAL1® NIR/ICG camera head should make sure that it is used in such an environment.			
Immunity tests	EN/IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) acc. to IEC 61000-4-2	± 8 kV Contact discharge ± 15 kV Air discharge	± 8 kV Contact discharge ± 15 kV Air discharge	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Rapid transient electrical Interferences/bursts acc. to IEC 61000-4-4	± 2 kV for power lines ± 1 kV for input and output lines 100 kHz repetition	± 2 kV/1 kV ± 1 kV for input and output lines 100 kHz repetition	Mains power quality should be that of a typical commercial or hospital environment.
Surges acc. to IEC 61000-4-5	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and fluctuations of the supply voltage acc. to IEC 61000-4-11	<u>Voltage Dip:</u> Dip to 0% for 1 cycle @ 0° phase angle Dip to 70% for 25/30 cycles @ 0° phase angle Dropout to 0% for 0.5 cycle @ 0°, 45°, 90°, 135°, 180°, 225°, 270° & 315° phase angles <u>Voltage interruption:</u> 100% for 250/300 cycles	<u>Voltage Dip:</u> Dip to 0% for 1 cycle @ 0° phase angle Dip to 70% for 25/30 cycles @ 0° phase angle Dropout to 0% for 0.5 cycle @ 0°, 45°, 90°, 135°, 180°, 225°, 270° & 315° phase angles <u>Voltage interruption:</u> 100% for 250/300 cycles	Mains power quality should be that of a typical commercial or hospital environment. If the user of the IMAGE1 S™ 4U Rubina, OPAL1® NIR/ICG camera head require continued operation during power mains interruptions, it is recommended that the IMAGE1 S™ 4U Rubina, OPAL1® NIR/ICG camera head be powered from an uninterruptible power supply or a battery.
Magnetic field at the supply frequency (50/60 Hz) acc. to IEC 61000-4-8	30 A/m at 50 Hz/60 Hz	30 A/m at 50 Hz/60 Hz	If image distortion occurs, it may be necessary to position the IMAGE1 S™ 4U Rubina, OPAL1® NIR/ICG camera head further from sources of power frequency magnetic fields or to install magnetic shielding. The power frequency magnetic field should be measured in the intended installation location to assure that it is sufficiently low.
Radiated, radio-frequency, electromagnetic field immunity test acc. to IEC 61000-4-3	3 V/m 80 MHz to 2.7 GHz *Refer table 2 for wireless Proximity RF field test levels	3 V/m 80 MHz to 2.7 GHz	
Immunity to conducted disturbances, induced by radio-frequency fields acc. to IEC 61000-4-6	3 V <sub>rms</sub> on 150 kHz to 80 MHz 1 kHz 80% AM modulation 6 V <sub>rms</sub> in ISM bands	3 V <sub>rms</sub> on 150 kHz to 80 MHz 1 kHz 80% AM modulation 6 V <sub>rms</sub> in ISM bands	

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**Table 3 – Test levels for Radiated and conducted Immunity Tests**  
Guidance and manufacturer's declaration – electromagnetic immunity

The IMAGE1 S™ 4U Rubina, OPAL1® NIR/ICG camera head is intended for use in the electromagnetic environment specified below. The user of the IMAGE1 S™ 4U Rubina, OPAL1® NIR/ICG camera head should make sure that it is used in such an environment.

Immunity tests	EN/IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Conducted RF disturbances acc. to IEC 61000-4-6	3 V <sub>rms</sub> 150 kHz to 80 MHz	3 V <sub>rms</sub>	Portable and mobile RF communications equipment should be used no closer to any part of the IMAGE1 S™ 4U Rubina, OPAL1® NIR/ICG camera head, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance: $d = 1.2 \sqrt{P}$
Radiated RF disturbances acc. to IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	where P is the nominal power of the transmitter in watts [W] according to the information provided by the manufacturer of the transmitter and d is the recommended separation distance in meters [m]. Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey <sup>a</sup> should be less than the compliance level in each frequency range <sup>b</sup> . $d = 1.2 \sqrt{P}$ 80 MHz to 800 MHz $d = 2.3 \sqrt{P}$ 800 MHz to 2.5 GHz  Interference may occur in the vicinity of equipment marked with the following symbol: 

Note: At 80 MHz and 800 MHz, the higher frequency range applies.

Note: These guidelines may not apply in all situations. The propagation of electromagnetic quantities is affected by absorptions and reflections of buildings, objects, and persons.

**a** Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the 4U Rubina, OPAL1® NIR/ICG camera head is used exceeds the applicable RF compliance level above, the 4U Rubina, OPAL1® NIR/ICG camera head should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the 4U Rubina, OPAL1® NIR/ICG camera head.

**b** Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

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<b>Table 4 – Emission class and group</b> <b>Guidelines and manufacturer's declaration – electromagnetic emissions</b>		
The IMAGE1 S™ 4U Rubina, OPAL1® NIR/ICG camera head is intended for use in the electromagnetic environment specified below. The customer or user of the IMAGE1 S™ 4U Rubina, OPAL1® NIR/ICG camera head should ensure that it is used in such an environment.		
RF emissions as per CISPR 11	Group 1	The IMAGE1 S™ 4U Rubina, OPAL1® NIR/ICG camera head uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause interference in nearby electronic equipment.
RF Emissions as per CISPR 11	Class A	The IMAGE1 S™ 4U Rubina, OPAL1® NIR/ICG camera head is suitable for use in all establishments other than domestic and those directly connected to the public low voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions as per IEC 61000-3-2	Class A	
Voltage fluctuations/flicker Emissions as per IEC 61000-3-3	Complies	

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<b>Table 5</b> <b>Recommended separation distances between portable and mobile RF communications equipment and the IMAGE1 S™ 4U Rubina, OPAL1® NIR/ICG camera head</b>			
<p>The IMAGE1 S™ 4U Rubina, OPAL1® NIR/ICG camera head is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or user of the IMAGE1 S™ 4U Rubina, OPAL1® NIR/ICG camera head can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the IMAGE1 S™ 4U Rubina, OPAL1® NIR/ICG camera head as recommended below, according to the maximum output power of the communications equipment.</p>			
<b>Rated maximum output power of transmitter [W]</b>	<b>Separation distance d [m] as a function of the transmitter frequency</b>		
	<b>150 kHz to 80 MHz</b> <i>d = 1.2 √P</i>	<b>80 MHz to 800 MHz</b> <i>d = 1.2 √P</i>	<b>800 MHz to 2.5 GHz</b> <i>d = 2.3 √P</i>
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23
<p>For transmitters rated at a maximum output power not listed above, the recommended separation distance D in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.</p> <p>Note: At 80 MHz and 800 MHz, the separation distance for higher frequency range applies.</p> <p>Note: These guidelines may not apply in all situations. The propagation of electromagnetic quantities is affected by absorptions and reflections of structures, objects and people.</p>			

**HF surgical instrument compatibility**

This equipment has been evaluated for compatibility with high-frequency surgical equipment. It has been tested according to IEC 60601-2-2 Annex BB.

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